

# **RAGAS DENTAL COLLEGE & HOSPITAL**

2/102, East Coast Road, Uthandi, Chennai - 600 119.

(A Unit of Jaya Group of Institution)



**FIRST B.D.S.**

**DENTAL MATERIALS PRACTICAL RECORD**

NAME : NAFIYA SABREEN . N.

# Ragas Dental College & Hospital

2/102, East Coast Road, Uthandi, Chennai - 600 119.

EXAMINATION REGN. NO. 541911551

DEPARTMENT OF DENTAL MATERIALS

Certified that this is the bonafide record of the work done

by Mr. / Miss NAFIYA SARREEN. N.

of First B.D.S. class in this Dental Materials Laboratory

at Ragas Dental College, Chennai.

During the academic year 20 - 20

*G. Kannan*  
14/12/20

Signature of the Staff

*[Signature]*

Signature of the  
Head of the Department

Submitted for the practical examination held in October / April 20

At Ragas Dental College, Chennai.

Date :

*[Signature]*  
11/12/20

*[Signature]*  
11/12/20  
Examiners

## GENERAL INSTRUCTIONS

1. Students must neatly dressed and wear their coats before entering the Laboratory.
2. Girls must have their hair braided and put inside the coat or have it tied-up to prevent mishap while working on Dental Motor.
3. Avoid wearing gold ornaments like rings and bangles when working with materials particularly with mercury.
4. Record books should be brought to the practical class without fail.
5. Necessary instruments should be brought to the practical class by every student.
6. Students should read the experimental procedure before coming to the practical class, perform the practicals carefully and record the observations.
7. Working place and the instruments should be kept clean.
8. All instruments and equipments should be carefully handled. Any damage to the laboratory equipments must be reported to the staff incharge.
9. Materials should be used sparingly during practical work.
10. The lids of the bottles and tubes should be replaced immediately after dispensing the materials.
11. Burners should be closed, and motors put off when not needed.
12. At the end of every exercise, the signature of the staff member should be obtained.
13. Record books should be preserved and submitted to the head of the department at the end of the academic session as they carry Internal Assessment marks.
14. All instruments can be procured through college Stores.

Sl No.	LIST OF INSTRUMENTS	Qty
1	BP Handle	1 No
2	Vulcanite Trimmer(Indian)	1 No
3	Carbide Burs for SHP 701, 702, 703	1 Each
4	Acrylic Trimmers	1 No
5	Sandpaper Mandril	1 No
6	Glass Slab	1 No
7	Cotton Holder	1 No
8	Chip Syringes with Nozzle	1 No
9	Glass Slide	1 No
10	Impression Trays (Reg)	4 Sets
11	Impression Trays (Reg)	4 Sets
12	Impression Trays (Reg)	4 Sets
13	Copying Pencil	1 No
14	Lecrons Carves(5)	1 No
15	Metal Scale	1 No
16	Bees and Paraffin Wax	1 Each
17	Blow Torch	1 No
18	Rubber Bowl S	1 No
19	Rubber Bowl L	1 No
20	Waxknife	1 No
21	Lecrons Carver	1 No
22	Wax Spatula	1 No
23	Plaster Spatula Straight	1 No
24	Plaster Spatula Curved	1 No
25	Cement Carrier	1 No
26	Bunsen Burner W/stopcock	1 No
27	Spirit Lamp	1 No
28	Flask No. 7	2 Nos
29	Handfile	3 Nos
30	Articulator Brown	1 No
31	Scissors Small Straight	1 No
32	Scissors Small Curved	1 No
33	Paint Brush	1 No
34	Glass Making Pencil	1 No
35	Metal Scale	1 No
36	Universal Plier	1 No
37	Glass Plate	1 No
38	Enamel Tray	1 No
39	MacIntosh Sheet 1 Metre	1 No
40	Hot Plate	1 No
41	Mouth Mirror With Handle	1 No
42	Straight Probe	1 No
43	Explorer/ Probe	1 No
44	Dental Tweezer	1 No
45	Plaster Knife	1 No

# CONTENTS

EX. NO.	DATE	EXPERIMENT	PAGE NO.
<u>I</u>	24/4/21	<b>Gypsum Products</b>	1
(i)	24/4/21	plaster of paris	6
(ii)	26/4/21	Dental Stone	8
<u>II</u>	15/03/21	<b>IMPRESSION MATERIALS</b>	10
(i)	28/4/21	Impression Compound	11
(ii)	30/4/21	Zinc Oxide eugenol paste	14
(iii)	30/4/21	Alginate	17
(iv)	3/5/21	Agar	22
(v)	8/5/21	poly sulphides	25
(vi)	8/5/21	Condensation Silicones	25
(vii)	4/5/21	Addition Silicones	31
(viii)	4/5/21	polyether	34
<u>III</u>	11/5/21	<b>DENTURE BASE RESIN</b>	37
(i)	12/5/21	Heat activated denture base resin	38
(ii)	12/5/21	Chemically activated denture base resin	43
(iii)	13/5/21	Light activated denture base resin	46
(iv)	13/5/21	microwave activated denture base resin	48
<u>IV</u>	27/6/21	<b>DENTAL CEMENTS</b>	49
(i)	28/6/21	Zinc phosphate Cement	51
(ii)	28/6/21	Zinc poly carboxylate Cement	53
(iii)	28/6/21	Zinc Oxide eugenol Cement	55
(iv)	28/6/21	Glass ionomer Cement	58
(v)	20/7/21	<b>DENTAL AMALGAM.</b>	61
(vi)	24/7/21	<b>DENTAL CERAMICS.</b>	65
(vii)	4/8/21	Composite resin.	67


## Gypsum Products

### Introduction:

- Gypsum was found in mine around the city of Paris. So, it is called as plaster of Paris. Mineral gypsum  $\text{CaSO}_4 \cdot 2\text{H}_2\text{O}$  is usually white in colour and it is formed as compact mass of gypsum. It is also an individual product.


### Manufacture:

The process of heating gypsum for manufacturing plaster is known as calcination. The mined gypsum is grounded and heated to a temperature of  $110^\circ\text{C} - 130^\circ\text{C}$  its loses off part of water by crystallization. It changes to  $\text{CaSO}_4 \cdot \frac{1}{2}\text{H}_2\text{O}$  remaining water of crystallization, is removed and the product formed as hexagonal anhydrite and orthorhombic hydrate.



### Setting Reaction

When plaster is removed, mixed in water, it takes up one and half molecule of water and becomes  $\text{CaSO}_4 \cdot 2\text{H}_2\text{O}$ .



This reaction is an exothermic reaction is same for all the gypsum products.

Three theories have been proposed to explain setting reaction, they are :-

- Colloidal theory
- Hydration theory
- Dissolution precipitation theory.

### Setting Expansion

The setting expansion of gypsum product can be detected during the change from hemihydrate to dihydrate during setting the gypsum product show linear expansion.

The setting expansion is due to crystallization mechanism, here, the enlargement of the dihydrate crystal which interlock with each other and obstruct the growth of a adjacent crystal.

If this process is repeated by thousand of crystal during growth, an expansion of entire

max. will occurs.

If the expansion is done underwater immersion it is called hygroscopic expansion of the setting if the expansion is without immersion under water is called normal setting expansion.

Types of gypsum products:

Type I - Impression plaster

Type II - Dental plaster

Type III - Dental stone

Type IV - Die stone

Type V - Dental stone, high strength, high expansion.

Type I - Impression plaster.

Impression plaster are composed of plaster of paris, where modifiers are added to regulate the setting time and setting expansion.

W/p water - 0.50 to 0.75

Setting time -  $4 \pm 1$  min

Uses :- It is restricted to final work impression of complete denture.

Type II - Dental plaster.

It is a type II plaster

W/p water - 0.45 - 0.50

Setting time -  $12 \pm 4$  min.



Uses :- It is used to fill the flask during flasking procedure in denture construction.

### Type III - Dental Stone

Intended for construction of cast in the fabrication of full denture that sets soft tissue. Stone dies are reproduction of prepared teeth as its compressive strength of 20-1 to 70 but does not exceed 24.1 mpa.

W/P ratio - 0.28 - 0.30

Setting time - 12  $\pm$  4 min

Use :- Dental Stone are preferred for cast and are used to process denture because the stone has adequate strength for that purpose and the denture is easy to remove after.

### Type IV - Die Stone

This material has high strength, hardness, abrasion resistance and minimum setting expansion.

A hard surface and sharp instrument are necessary for a die stone.

W/P ratio - 0.22 - 0.24

Setting time - 12  $\pm$  4 min

Use :

A hard surface is necessary for a die stone because cavity preparation is filled with wax and carried flush with margins of die. A sharp instrument is used, stone must be resistant to abrasion.

Type II - Dental stone, high strength, high expansion.


This material has higher Compressive Strength. This is obtained by lowering w/p ratio.

w/p ratio = 0 - 18

Setting time =  $12 \pm 14$  min

Use :

The rationale of increase in setting expansion limit is that certain newer alloys such as base metal, have greater casting shrinkage than the traditional noble metal alloys. So, high expansion is required for fabrication of die to aid in compensating for the alloy solidification shrinkage.



## Plaster of Paris

### Introduction



Model plaster is used to fill the flask in denture construction, it reproduces minute detail, while pouring the cast an impression. This is very hard and strong.

### Composition

- $\text{CaSO}_4 \cdot \frac{1}{2} \text{H}_2\text{O}$
- $\text{CaSO}_4 \cdot 2 \text{H}_2\text{O}$
- Modifiers (accelerators & Retarders)
- Common accelerator is  $\text{K}_2\text{SO}_4$ .
- Retarder - Borax
- Surface active agent

### Manufacturing of Plaster of Paris:

Gypsum is ground in a open container and heated at  $110^\circ\text{C} - 130^\circ\text{C}$  by a process called dry calcination. Fine crystals aggregate of fibre to form a B-type as B-hemihydrate forming an end formed product of needle like crystal.



### Manipulation:

Instruments required: Flexible rubber bowl, Soft powder mixing Spatula, Proper w/p ratio is assumed.

### Procedure :-

- Measure quantity of water in clean and dry rubber bowl is taken.
- Pre-weighed quantity of powder is shifted into the bowl.
- It is allowed to settle for 30 Sec to avoid air entrapment.
- Spatulation is done with a hard Spatula vigorously for 45-60 Sec till smooth and Creamy mix is obtained.

w/p Ratio : 0.45 - 0.50

Working Time : 3 min

Setting Time :  $12 \pm 4$  min

### Properties :-

- It is in a white natural colour.
- Relatively weak as Compressive Strength - 9 mpa and tensile strength - 0.6 mpa
- Low w/p ratio for high strength.
- Dimensionally accurate and stable.
- Good colour contrast.

Uses :-

- Used for making diagnostic Cast and model
- for mounting Cast on articulator


## Dental Stone

### Introduction:

Also called as type III gypsum product used for reproduction of die of prepared teeth for fabrication of prosthesis and for construction on full dentures.

- $\text{CaSO}_4 \cdot \frac{1}{2} \text{H}_2\text{O}$  - Principle Constituent
- Surface active agent
- $\text{K}_2\text{S}_2\text{O}_8$  - accelerator
- Borax - retarder

Gypsum is calcinated under steam pressure in a autoclave or a kettle at  $120^\circ\text{C} - 130^\circ\text{C}$  for 5-7 hrs. The product obtained is much stronger and harder than  $\beta$ -hemihydrate called  $\alpha$ -hemihydrate.



Crystals in the form of rods and prisms is seen microscopically.

### Manipulations:-

This involves proportional mixing and a dry gypsum product and water to get a Smooth Creamy mix of Stone.

### Materials required

Para bolt flexible rubber band, stiff spatula

W/p ratio : 0.85 - 0.90

Working time : 3 min

Setting time : 12 + 14 min

### Properties:-

- It is green or yellow in colour.
- Comparatively smaller in size.
- Compressive strength - 20.7 mpa - 34.4 mpa
- Setting expansion - 0.6% - 0.8%
- Hardness - 80
- Easy to manipulate
- Cheap and economic
- Prothemic in nature
- Absorb moisture when placed in tissue.

### Uses:-

It is used to construct a cast for the fabrication of complete denture. It is used for construct study models and for preparing master cast and models.

✓

## Impression Materials

### Introduction:

Impression material are used to register or reproduce the form and relationship of teeth and oral tissue in a negative aspect.

The negative reproduction of tissue is called Impression.

### Classification:

classified based on mode of setting and elasticity (or) plasticity of material.

By chemical/irreversible

By physical/reversible

Elastic

Alginate / Hydrocolloid  
Elastomeric material

Agar hydrocolloid

- Polysulfide
- Condensation silicone
- Addition silicone
- Polyether

Rigid

Impression Plaster  
Zinc oxide eugenol

Impression Compound  
Impression wax



# Impression Compound

## Introduction:

Impression Compound is one of the class of the impression materials.

ADA Specification number: 3. It is also called modelling plastic mainly used for making impression of edentulous ridges.

## Classification:

- Type I - Impression Compound
- Type II - Tray Compound

## Mode of Supply :-

Supplied as sheets, sticks, cakes and cones in various colours.

## Composition :-

- Resin (30%) - Major ingredient
- Copal resin (30%) - Provides thermo plasticizer.
- Carnauba (10%) - provides thermo plasticizer.
- Acid (5%) - plasticizer.
- Talc (25%) - filler.
- Colouring agent.

## Manipulation :-

The compound is placed in hot water, held for minute and it is taken out. It is manipulated by wet kneading method.



### Materials Required :-

Hot water bath - impression Compound.

Two Method - Dry heat and hot water bath

**Dry heat** - Small amount of Compound are obtained over a flame directly.

**Hot water bath** - The Compound is immersed in hot water bath of  $65-70^{\circ}\text{C}$  and above and softened into a soft dough kept on tray and impression is taken and set to harden, withdraw from mouth.

### Properties :

- The fusion temperature of the - impression Compound is  $42.5^{\circ}\text{C}$ .
- The thermal conductivities are very low.
- The Co-efficient of thermal expansion -  $0.3\%$ .
- The dimensional stability of impression Compound is poor. Surface detail product is relatively less because of its high viscosity and low flow.
- It is non-toxic and non-irritant.
- slow flow is needed during impression making.

### Advantages :-

- The Materials can be reused for a number of times.
- It is cheap and Economic.

- The Accuracy can be improved by flaming the surface and moulding it with oral tissue.

### Disadvantages :-

- Compress the soft tissue.
- Distortion due to less dimensional stability.
- Difficult to remove if there is severe undercut.

### Uses :-

- To make an impression of complete edentulous
- It is used in border moulding technique for fabrication of complete denture.
- To make impression of individual teeth.

~~Handwritten mark~~

## Zinc Oxide eugenol impression paste

### Introduction:-

Zinc Oxide eugenol is properly known as the impression material for making impression of edentulous arches for the construction of complete denture.

### Classification:-

Type I - Hard

Type II - Soft

### Available as:-

- Base plate - white colour.
- Accelerator plate - red colour.

### Composition:-

#### Base plate:-


- Zinc Oxide - 37%
- Vegetable Oil - 13%

#### Accelerator paste:-

- Clove Oil (or) eugenol - 12%
- Green (or) polymerised resin - 13%
- Filler - 20%
- Resin balsam - 10%
- Calcium chloride of colour - 5%

### Setting Reaction:-

The setting reaction is typical acid base reaction to form chelate. The reaction is known as chelation and product is Zinc Sulfate.



### Manipulation:

The mixing is done on a non-sticky oil paper or glass slab. Two scooper of paste of same length and width are from each tubes are squeezed into the mixing slab and a flexible stainless steel spatula are mixing until the uniform colour is obtained.

mixing time :- 1 minute

Setting time :- Initial Final

Type I 3-6 min 10 min

Type II 3-6 min 15 min

### Properties:

Consistency and flow: - good flow and homogeneous consistency

Type I - 30-50 min

Type II - 30-45 min

Compressive Strength - 7 mpa.

### Advantages :-

- Register accurate surface detail.
- It is dimensionally stable.
- Minor defect can be corrected locally without discarding a good impression.
- Does not require separating medium.

### Disadvantages :-

- It is sticky in nature.
- It requires special tray for impression making.
- Eugenol causes tissue irritation and burning.
- Not used for making impression of teeth and undercuts as it is elastic.

### Uses :-

- Used as final impression material in completely edentulous patient with few undercuts.
- Temporary cement are used to relieve an area to improve accuracy of record Surgical pattern.

### Disinfection :-

Zinc Oxide eugenol can be disinfected by immersion in 2% glutaraldehyde in 1:2:3.  
The manufacturer recommendation for proper disinfection should be followed.

✓

## Alginate

### Introduction :-

Alginate is inconvertible hydrocolloid. It is process of cross linkage algin acids with Calcium ions.

ADA Specification - 18

Type I - Fast Setting

Type II - Normal Setting.

### Composition :-


Supplied in Container as paste form.

Ingredient	Percentage	Uses
Na alk alginate	15	Main ingredient
Calcium Sulfate	16	Reaction
Zinc Oxide	4	Fillers
Potassium titanoaluminate	3	Cuprous hardener
Diatomaceous earth	60	filler
Calcium phosphate	2	Retarder
Colouring and odouring agent	Trace	Flavour addition

### Setting reaction :-

The typical sol gel reaction can be described simply as a reaction of soluble alginate with calcium sulfate and formation of an insoluble calcium alginate gel. Cases react rapidly to produce insoluble calcium alginate in an aqueous solution. The production of calcium alginate does not

allow significant sufficient working time. So, water solution salt such as sodium phosphate is added to solution to prolong working time, when the supply of trisodium phosphate is exhausted. The calcium ion begin to react with potassium alginate to produce calcium alginate as follow. The retarder must be adjusted carefully to provide to proper gelation time.



### Manipulation:-

The powder supplied in a can is scooped by Shaking and 1 Scoop of powder is used for the one measure of water with the measuring Component supplied along with pack.

Cooling the water increase the working time and setting time.

The powder is sifted into water and it is mixed with help of Curved Spatula in a flexible bowl by a figure of eight method.

W/P ratio : 15 g 40 ml

Mixing time : Fast Setting : 45 Sec

Normal Setting :  $\pm$  60 Sec

Setting time : Fast Setting : 75 Sec

Normal Setting : 120 Sec

### Properties Flexibility :-

It is about 14% and act as stress of 100 gm/cm<sup>2</sup> same hard set material have values from 5% to 8%, lower w/p ratio result in lower flexibility.

**Elasticity :-** Elasticity is less when compared with agar and about 97.3% elastic recovery. Some permanent deformation is less than when impression is removed from mouth quickly.

**Strength :-** Compressive Strength - 500 - 8000 g/cm<sup>2</sup>  
Tear Strength - 30 - 700 g/cm<sup>2</sup>

### Advantage :-

- Inexpensive
- Easy to manipulate
- Pleasant taste and Odour
- Hydrophilic Easily paired with dental Stone
- They can be used with Stock trays.

### Disadvantages :-

- Tears easily.
- Must be paired immediately after removal.
- Has limited detail reproduction.
- Can be used for single cast
- Incompatible with die material



### Disinfection:

Disinfecting the hydrocolloid impression by using household bleach, iodophores or synthetic bands.

It is used thoroughly the disinfection is sprayed liberally on surface. It is immediately wrapped in disinfectant-soaked paper towel & placed in a sealed plastic bag for 10 mins.

### Uses:

Routine impression for diagnostic and incomplete edentulous and dentulous patient.

1/24/21

## Agar

### Introduction:

Agar is an organic hydrophilic colloid prepared from certain types of seaweeds. It is a sulphonic ester of linear polymer of galactose. The principle ingredient by weight in water, greater than 80%. Modifiers are present in minor amount and exert a considerable influence on properties. ADA specifications number: 1

### Composition:

Ingredients	Percentage	Function.
Agar	13-17%	Main ingredient.
Borax	0.3 - 0.5%	Improves strength.
K <sub>2</sub> SO <sub>4</sub> retarder	1-2%	Retarder.
Thermoplastic material	0.5 - 1%	Plasticizer.
Alkyl Benzoate	0.3 - 0.5%	Plasticizer.
Colouring & flavouring	Trace	Colour & flavour.
Water	Balance	Dispersion medium.
Hard Wax	0.5 - 1%	Filler.

### Manipulation:

- Prior to impression procedure.
- Impression making.
- Post impression procedure.

### Preparation of material:

Special equipment for liquifying and rearing agar material is essential. The hydrocolloids are usually supplied in syringe and tray material. The only difference between them is the colour and viscosity of material. Tray material is more viscous.

Reversal of hydrocolloid from gel to sol state. A special cooker with 3 chambers and cooling trays are used.

Heat in water bath at  $100^{\circ}\text{C}$  for 8-12 min & stores in second chamber at  $46^{\circ}\text{C}$  water for 2-10 min. After heating the tray in oral cavity. Cool with water. No less than  $13-20^{\circ}\text{C}$  for gelation. After impression is removed from the mouth, wash to remove debris to avoid interference during setting of Gypsum.

After initial setting of gypsum, model & impression in a humid location.

### Properties:

**Flow:** The material is sufficiently fluid to record fine detail if it is correctly manipulated.

**Gelation temperature:** The agar solution should be reheated to get at a temp. of  $37^{\circ}\text{C}$  but not more than  $45^{\circ}\text{C}$ .

**Flexibility:** The flexibility of gel from time of removal is between 4-12% when a stress of  $4.212\text{ps}$

is applied.

Elastic Recovery: highly elastic, extend upto 98-97.

Cooking Time: 7-15 min

Setting Time: 5 min.

Advantages:

- Inexpensive
- Have pleasant colour
- Non toxic and non staining
- Doesn't require custom tray as adhesive
- Component do not require mixing.
- material is hydrophilic and can be used in presence of moisture.

Disadvantages:

- Material requires use of expensive equipment.
- Tear Easily
- Must be poured immediately
- Dimensionally unstable.

Use:

Used as duplicatory material in the construction of cast partially denture cast the accuracy for precisely cast and other restoration.

Disinfection: Agar hydrocolloid are disinfected by immersion in sodium hydrochloride, iodophor phenolic glutaraldehyde.

## Polyulfides

### Introduction:

Polyulfides are widely used for woven and bridge fabrication due to high strength and relatively unimpromen and unimpromes curing time is needed They are supplied in tubes of base plate and catalyst paste are mixed together. It is also called as thiol. ADA-19.

### Composition:

#### Base Plate:

Liquid polyulfate polymer - 80 to 85%

Trest filler - 16-4%

Plasticizer.

#### Accelerator Paste:

Lead Dioxide - 60-65%

Diethyl Ophthalate - 30-35%

Sulfate - 3%

Magnesium Retarder - 2%

### Manipulations:

The catalyst and base plate is dispersed in glass slab of same length. The catalyst paste is collected in stainless steel spatula.

The mixture of spread over the mixing pad. The mass is then trapped for 10th block and spread back

and forth on mixing pad. This process is continued until mixing pastes are uniform in colour.

### Setting Reactions:

Lead Dioxide react with polysulfide polymer & causes

- Chain lengthening by oxidation of terminal SH groups

- Cross linking by oxidation:



### Properties:

Unpleasant colour & colour. Excellent reproduction of surface details. Shelf life is good. The mixing shrinkage is high - 0.45%. It has high permanent deformation and its 35% average elasticity.

### Advantages:

- Long working time.
- Good tear strength
- Good flow
- Good production & surface detail.
- High flexibility.

### Disadvantages:

- Needed use of custom tray
- Bad colour.
- Lower viscosity

- Lead dioxide material stains clothing

Use:

Polysulphide Impression material is often used with custom tray to increase the accuracy of impression. Light and medium body material are used for complete denture impression.

Disinfection:

Disinfected by sodium hypochloride, sodium phosolite, glutaraldehyde

## Condensation Glucone.

### Introduction:

They are used mainly for crown and bridge impression making they are ideal for single unit. Inlay supplied by two pastes - as paste liquid catalyst system.

### Composition:

#### Base Paste:

Polydimethyl Siloxane - Principle constituent

Dibutyl Oxthalate - Plasticiser.

Colloidal silica - Filler.

Platinum - 33-75%.

#### Catalyst Paste:

Metal Organic Ester.

Tin Acetate.


Dibutyl tin Hexanoate.

### Manipulation:

Equal lengths of base & accelerator is extended into mixing pad side by side without touching.


The accelerator paste is incorporated into base paste. Just before shading the tray, the material should be spread as thin layer dimensional network is formed. Stannous acetate acts as catalyst.





Setting Reaction:

It is a condensation reaction occurs as a result of cross linking between orthoethyl silicate & terminal hydroxyl group of dimethyl silicone to form Si-O network.



Mechanical Properties:

Accuracy is greatly improved by first taking impression with a highly filled silicone putty and after setting taking a second impression with light bodied silicone. Thus the final total shrinkage is lower.

Advantage:

- Clean, pleasant material for patient.
- Highly elastic and setting time can be controlled with amount of accelerator.
- Eliminates need for custom tray.

### Disadvantage:

Brainwate due to shrinkage on standing and should be poured within 1hr. They are very hydrophobic, require a dry field and are difficult to pour in stone.

### Uses:

Heavy and light bodies syringe material may be used to make condensation silicone impression using the multiple mix technique described for poly sulphide.

### Disinfection:

Disinfectant by sodium hypochlorite, iodophor complex phenolic glutaraldehyde.

## Addition Silicone.

### Introduction:

Represent an advancement in curable resin condensation silicone. This has been achieved by change in polymerization reaction to an addition type and elimination of no alcohol by product that takes place.

### Composition:

#### Base paste:

Hybrid Silicone.

Palladium.

Silicone polymer OH group.

#### Accelerator paste:

Potassium salt.

Silicone polymer with vinyl group.

### Manipulation:

It is available only in two part system. The two systems are dispersed in volume by equal no. of scoops of each part. The best mixing technique is to mix with finger until a uniform colour is not obtained during manipulation that it is not uniform and it causes distortion.

### Setting Reaction:

It is an addition reaction. In this case the polymer is terminal with vinyl group and linked with siloxane. This reaction is activated by platinum salt.



### Properties:

- \* Setting time - 4-7 min.
- \* Indenting time - 45 sec.
- \* Good shear strength.
- \* Low flexibility and shelf life ranges from 1-2 years.

### Mechanical Properties:

The working and setting time of addition siloxane are as faster than polysulphide rubbers. It often used to stand the working time.

They have excellent elasticity and show very less dimensional change upon storage. Therefore can be safely poured later.

### Advantages:

- High strength and dimensional stability.
- Resisting from deformation.
- Pleasant colour.

### Disadvantages:


- Expensive, more rigid than condensation silicone and difficult to remove from undercut.
- Supports the latex, gloves & rubber can exhibit polymerization.

### Uses:

- Bite Registration material.
- Diagnostic aids.
- Fabricate mold for fabrication of temporary crown, heavy viscosity key material.

### Disinfection:

Immersed in sodium hypochlorite, iodophor, phenolic gluteraldehyde.



## Polyester.

### Introduction:

Used for accurate impression for a few prepared teeth without severe undercuts.

Due high stiffness and short working time restrict their use to impression of a few teeth.

### Composition:

#### Base Plate:

Polyester Polymer

Colloidal silica

Glycol ether (or) phthalate

Accelerator paste

Phosphoric sulphate ester

Colloidal silica filler

### Manipulation:

Equal length of base & catalyst paste is dispensed in groove of slab with graduation. The catalyst paste is collected in spatula and is distributed over the base. The mixture is then scrapped back and forth over the mixing pad.

### Setting Reaction:-

It is a reaction between benzoinic ring and an end of branched polyester molecule. The main chain is a copolymer of ethylene dioxide and tetrahydrofuran.

Cross linkage is brought about by automatic Sulphonic ester.

polyether + Sulphonic acid ester  $\rightarrow$  Cross linked rubber.

Properties :-

- Setting time - 5 min
- Mixing time - 30 Sec
- Hydrophilic So moisture in impression is not critical.

Mechanical Properties :-

Because this rubber absorbs water and changes dimension. Storage in water is not recommended. Elastic recovery values average 80-85% between base poly sulphides and addition.

Advantages :-

- Pleasant to handle, easy to mix.
- These materials are more accurate than polysulfide on condensation silicone impression.

Disadvantages :-

- Include high cost, short working and setting time, high stiffness (after setting which limit their use).
- Their bitter taste is objectionable to some patients.
- They cannot be kept for long time in disinfection solution.

Uses :-

Used in dispensable staple using strategies of  
offered partial denture, medium viscosity used in  
border moulding and thin viscosity as seal impressions.

Disinfection :-

By Immersed in Sodium hypochlorite.

~~Unit~~



## Classification of denture base resins

### Denture base material



#### Metallic

- Cobalt Chromium
- Gold alloys
- Aluminium Stainless Steel

#### Non-Metallic

- Styrylic Resin
- Vinyl Resin

### Denture base material



#### Temporary

- Self cure acrylic
- Shellac base plate
- Base plate wax

#### Permanent

- Heat cure wax
- Light cure wax
- Pain type wax
- Metallic wax

## Heat activated denture base resin

### Introduction :-

Heat activated material are used in the fabrication of nearly all denture base. Thermal energy required for polymerisation of such material should be provided using a waterbath or microwave Oven. Because of the prevalence of microwave Oven, herein emphasis is placed on heat activated system. ADA - 42.

### Composition

	Ingredients	Function
Liquid	Methyl methacrylate Di butyl phthalate Glycol methacrylate Hydroquinone	Plasticizer plasticizer prevent Inhibitor of setting
powder	Compound of mercury Sulfide Zn or TiO <sub>2</sub> Di butyl phthalate Dyed Organic filler Inorganic particle like glass fibre	Dyes  Plasticizer

### Manipulations: -

- preparation of waxed denture pattern.
- preparation of Split material
- application of Separating medium, mixing of powder and liquid of denture base resin.
- packing
- Lining
- Coding
- Deflasking, finishing and polishing.

### Preparation of waxed denture pattern :-

The teeth is arranged in wax (or) wax pattern is oflashed in the denture flask. The coat is placed in base of flask and it is invested with plaster of Paris. It is called as base oflasking.

### Preparation of Split mould :-

After the plasterizing under a second mix of plaster is replaced, prepared and poured into the denture flask. The investment is added until all Surfaces of tooth arrangement and base is Covered Completely. The Remainder of all the flask is filled with lid is gently placed in position. The flask is placed in the removed from water and Segment are Separated.

### Application of Separating Medium :-

Currently most popular Separating medium agent are water Soluble alginate Solution. When Applied to dental Stone Surfaces leads dissolution produce them.

relatively insoluble calcium and alginate films that prevent direct contact of denture base resin and surrounding dental stone.

Mixing of powder and liquid :-

Polymer of monomer ratio 3:1 by volume and 2:1 by weight. The measured liquid is poured into clean, dry, mixing bowl, powder is added and allowing each of particle is wetted by monomer.

Packing :-

The packing process should be performed while the denture base resin is dough stage. The resin is removed from its mixing container and roller it in a mixing container and roller it in pipe form (a) rod from which is bent into horse-shoe shaped and placed in portion of flask that house prosthesis teeth.

Polymerisation Cycle :-

proceed at 74°C approx 2 hours and then increased to 100°C for 1 hr. The flask is cooled after the procedure for 30 min then immersed in cold water for 15 min.

Setting Reaction :-

Setting Reaction takes place between polymer and monomer is achieved by heat and pressure.



### Wet Sand Stage :-

polymer gradually settles, into monomer forming a fluid or coherent mass.

### Sticky Stage :-

Monomer attack the polymer, it become penetration into polymer. The mass is sticky and strong when touched and pulled against.

### Dough Stage :-

As monomer diffuse into the polymer it become more smooth and dough like and it doesn't adhere to walls of container an idea stages of making DGE is a demand must to rise.

### Rubbery Stage :-

Monomer disappear by penetrating into polymer. Non elastic and can be molded at this stages.

### Ideal Properties :-

Good Appearance, high flexural strength, high impact strength, high stiffness, long fatigue life, high creep resistance, high creep resistance, high rigidity.

### Advantages :-

- High glass transition temperature.
- Easy of fabrication.
- Low capital cost.
- Good surface finish • Good Appearance.

### Dis Advantages :-

• Free monomer, Content of formaldehyde can cause sensitization, low impact strength and radiopacity. Flexural strength low enough to penalize poor denture design, fatigue life is short.

### Uses :-

- Fabrication of permanent denture base in implant supported over denture.
- Used in occlusal splints.

## Chemically activated denture base resin

### Introduction:-

In Addition to the use of heat energy Chemical activation are also used in indirect denture base polymerization. It does not require application of heat and this is completed at room temperature.

### Manipulation:-

- Manipulated by following technique.
- Sprinkle on technique.
- Adapting technique.
- Fluid Resin technique.
- Compression moulding technique.
- Injection moulding technique.

### Sprinkle on technique:-

Operating medium is applied on cast and powder and liquid is alternatively from dropper. The prosthesis is build section by section until the completion to improve the strength and appearance is used in hot pressure to laminate it.

**Adapting technique:-** Powder and liquid is proportionally mixed in porcelain jar. It under dough stage. It is quickly removed and adapted on cast and manually moulded.

### Fluid resin technique :-

A completed arrangement of teeth positioned in a flask and is inserted with agar. The flask is removed from water bath and sprues are attached for resin introduction. The resin is poured through sprue and completely protrusion is removed.

### Compression mould technique :-

Completed tooth arrangement is prepared for flasking procedure, the waxed cast is embedded in freshly poured stone. The surfaces of tooth are covered with stone.

### Injection mould technique :-

Completed tooth arrangement is prepared from flasking procedure, the sprues are placed for introduction of resin, wax elimination is done, no separating medium is applied on mould space, and resin is introduced into the cavity and placed in water bath.

### Setting reaction :-



To obtain a workable mass it has to undergo 5 stages. In Steady Stage the mixture is coarse or grainy. In Stringy Stage, the monomer enters bead stage increased polymer chain. enters the solution and dissolves and enters a clough Stage.



### Advantages :-

- Easy to deflex.
- Dimensional accuracy.
- Capable of higher flexural strength.

### Disadvantages :-

- Not cheap.
- Increased creep.
- Colour instability.
- Reduced stiffness.
- Tooth adhesion failure.

### Uses :-

- Used for relieving and repair work.
- Used in diagnostic aid.
- In the construction of immediate denture.
- In maxillofacial prosthesis.

# Light activated denture base resin

## Introduction:-

Light activated denture will compare well to conventional heat and materials and some of impact strength and hardness, but low elastic modulus.

## Composition:-

Ingredients	Uses
Urethane dimethacrylate	Principle ingredient
microfine silica	filler
Camphor quinone	Photo initiator.

## Manipulation:-

It is polymerised in a light chamber with UV to blue light of 400-500 nm, an high intensity quartz halogen bulb. The denture is rotated continuously to provide uniform exposure to light.

## Setting reaction:-

Light a activator. The Stage of polymerisation reaction are wet, Sandy, Sticky, dough, rubber and stiff stage.

### Advantage :-

The finished product is durable and light weight. Speed of polymerisation is accomplished.

### Dis Advantages :-

It require special instrument for manipulation temporary denture base.

### Uses :-

- Used for denture base.
- for making special tray.
- Used for making obturator.
- formation of light activated denture base Contact Copolymer methac diacrylate and acryle along with silica filler.

Microwave activated denture base resin

poly(methyl methacrylate) may be polymerised using microwave energy. This technique employs specially formulated resin in non-metallic flask.

4/2

## Dental Cement

### Introduction :-

Dental Cement are material of different use including restoration cutting and thin section. They are generally material of comparatively low strength but have extensive use in dentistry. The first cement of dentistry was introduced in 1755 by Singer who created Zinc Oxide cement. Sorel

### Classification :-

I) Iso - Standard Covering Cement.

Iso 9917-1: 2003 - water based cement part-I  
Powder liquid base cement.

Iso 9917-2: 1998 - water based cement part-II  
liquid activated cement.

Iso 9909 - 2004 - Zinc oxide eugenol and  
Zinc Oxide non eugenol cement.

Iso - 4029 - 2000 - polymer based filling,  
restorative and luting  
material.

ii) According to setting reactions:-

- Acid-base reaction cement
- polymerising cement
- Dual cure cement
- Tissue cement

Ex) Classification of Cement based on application:

Type I :- Luting Agent. Eg: Zinc phosphate, fine grain for cementation, luting medium grain for base, Orthodontic appliance.

Type II :- Restorative application. Eg: Glass ionomer Cement

Type III :- Lining or Base application  
Eg: Calcium hydroxide.

General Structure :-

On mixing powder and liquid only a part of powder react with liquid and form a dental material with a

- a core of untreated powder surrounded by the reacted material.
- a matrix of reacted material.

Uses :-

- Final Cementation

- Temporary Cementation

- Bases

- Long term restoration

- Temporary and intermediate restorations

- Pulp therapy

- Obturator

- Liner

- Root Canal Sealants.

## Zinc phosphate Cement

### Introduction :-

Zinc phosphate is the oldest of cutting agent for the term Crown and bridges. Also used for Amputation.

### Applications :-

- Luting of restoration
- High strength base
- Temporary restoration
- Luting of Orthodontic band and brackets.

### Classification :-

Type I - Fine grained for luting.

Type II - Medium grained for luting.

### Composition :-

Powder :- Zinc Oxide (90.9%) - Principle Constituent  
 Magnesium Oxide (8.2%) - Aide in Sintering  
 Alumina Oxide (0.2%) - Improves Smoothness  
 Silica (1.4%) - Filler aide in Sintering.

### Liquid :-

phosphoric acid (35.2%) - Reacts with Zinc  
 water (36%) - Control the rate of reaction  
 Aluminium phosphate (10.5%) - Buffer with Zppos.  
 Zinc (4.1%)  
 Aluminium (2.2%)

## Setting Reaction:-

When powder is mixed with liquid phosphate Zn ions. The aluminium particles and release of cementation like aluminium complex with phosphate acid Zn ions to form Zn aluminium phosphate gel.

## Manipulation:-

Powder is taken on slab along with liquid. Liquid is added to powder in small increments, mixing is done in brush circular motion.

Mixing time - 1 - 1.5 min

## Properties:-

- Compressive Strength - 100 - 119 MPa
- Tensile Strength - Set cannot be little
- Modular of Elasticity - 13.7 GPa
- Solubility -  $10^{-10}$  Solubility
- Biological property - pH acidic pulp response is moderate

Advantage:- Long track record

- Good Compressive strength.

## Disadvantages:-

- No chemical adhesion
- No anti-cariogenic effect.
- Pulp irritation
- Poor esthetics.



## Zinc Polycarbonate Cement

### Introduction :-

Canadian biochemist Smith developed the first polycarbonate cement in 1968 by substituting the phosphoric acid of zinc phosphate cement with polyacrylic acid.

### Application :-

- For luting permanent restoration.
- As base for liners.
- Used in orthodontics for cementation of band.
- Used in root canal filling in endodontics.

### Available as

- Powder and liquid
- Water Settable Cement
- As pre-capsulated powder and liquid.

### Composition :-

powder :- Zinc Oxide - Base Ingredient  
Magnesium - Principle modifying Oxides  
of barium and aluminium - Small amount.

Stannous fluoride - Anti-Cariogenic

### Liquid :-

Aqueous solution of polyacrylic acid or polymer of acrylic acid with unsaturated carboxylic acids.

### Manipulation :-

powder & liquid are taken in cooled glass slab. powder is incorporated in bulk into liquid. mixing is done using cement spatula.  
Mixing time - 30 - 40 Sec.

### Setting Reaction :-

when powder and liquid are mixed, the surface of powder particles are attacked by acid releasing zinc, magnesium, tin ions. These bind to polymer chain in Carboxyl group. They also react with Carboxyl group of adjacent polyacid chain to form cross linked to balls.

Setting time :- 7-9 min.

### properties :-

- Compressive strength - 55 mpa
- Tensile strength - 6.2 mpa
- Solubility - absorb water and slightly soluble.
- Biocompatibility - pH - 1.0 - 1.5.
- pulp response is mild
- Anti - Caries

### Advantage :-

Less irritation to pulp.  
Chemically bonds to tooth.

### Disadvantage :-

Limited to fluoride release.

## Zinc Oxide eugenol cement

### Introduction:-

Cement have been used extensively in dentistry since the 1870's. The general 9 Cement are of low strength. They are least irritating and known to have antiseptic effect.

### Classification:-

- Type I - Temporary restoration.
- Type II - permanent restoration.
- Type III - Bases and temporary restoration.
- Type IV - Caring jaws.

### Available as:-

- powder and liquid
- Two paste systems

### Compositions:-

#### Powder:-

- Zinc oxide - 69% - principle ingredient
- white zinc - 29.3% - Reduce brittleness
- Zinc stearate - 1% - Accelerator plasticizer.
- Zinc acetate - 0.7% - Accelerator improve strength, Magnesium oxide.

#### Liquid:-

- Eugenol - 25% - react with ZnO
- Oil of - 15% - plasticizer.

Manipulation:

Measured quantity of powder and liquid is dispersed into cool glass slab. The bulk of powder is incorporated with a liquid spatulated thoroughly in circular motion with a stiff backed Stainless Steel Spatula. A bit of orange is used to clean enamel cement from instrument.

Setting Reaction:

Fast hydrolysis of Zinc Oxide takes place and form zinc hydroxide later this react with enamel to form zinc enamel.



Setting time - 4 - 10 minutes.

Properties:


- Compressive Strength - 5 - 5.5 mpa
- Tensile strength - 0.32 - 5.32 mpa
- Thermal property = Excellent
- High Solubility
- pulp response - mild.

## Modified Zinc Oxide eugenol Cement

Aluminic modified cement. It has greater strength and used as intermediate, filling material and as base. Liquid is substituted with orthophosphoric acid. Alumina is added to powder, polymer reinforced Zinc Oxide eugenol cement.

## Special Zinc Oxide eugenol product.

ZOE is very popular as endodontic sealer also traditional formulation along with gutta percha resin. These material are used in endodontics therapy to root canal.



## Glass Ionomer Cement

### Introduction:-

Glass ionomer Cement evolved from general dissatisfaction with silicate Cement. The first usable c/p.c. was formulated in 1972 by Wilson and was known as aluminosilicate.

### Application:-

- Anterior esthetic restorative material for class IV
- Restorative material for eroded areas and class III restoration
- As a luting Agent
- As a liner and base.
- For core build up
- To a limited extension on pit and fissure Sealant.

### Classification:-

- Type I :- for luting
- Type II :- for restoration
- Type III :- for liners and bases.

### Available as:-

- Powder / liquid in bottle
- preportioned powder / liquid in bottle.
- Light cure System.
- powder / distilled water (water settable)

Composition :-

powder - silica ( $\text{SiO}_2$ ) - 41.7%  
 Alumina ( $\text{Al}_2\text{O}_3$ ) - 45.6%  
 Aluminium fluoride - 1.6%  
 Sodium fluoride - 1.5%  
 Aluminium phosphate - 8.3%

Liquid :- 50% aqueous solution of polyacrylate acid in the form of co-polymer with stearic acid, maleic acid, iteicarbonylic acid tend to increase relatively of fluid decrease the viscosity and reduce tendency of gelation.

Manipulation :-

powder and liquid are dispersed onto a coat absorbant paper pad. The powder is added in two or more iteicarbonylic acid tend to increase relatively of fluid decrease the viscosity and reduce tendency of gelation. Spatula is used for mixing and following technique is used.

Setting reaction :-

- bonding when powder and liquid are mixed together, the acid attacks glass particles, thus calcium, aluminium, sodium and fluoride.
- Calcium cross link - The initial set occurs when calcium cross link with polyacrylic chain.

- Aluminium cross link with polyacrylic acid chain  
Sodium and fluoride ions - Sodium ion may replace hydrogen ions in - Carboxylic groups. The rest combine F its form NaF.
- Hydration initially water serves as medium, but it hydrates the matrix adds strength to the silica gel sheath.
- The unreacted glass particles is sheathed by silicic acid.

### Properties :-

- Compressive Strength - 85 - 105 MPa
- Tensile Strength - 62 - 66 MPa
- modulus of elasticity - 7.3 GPa
- Solubility - Initially high
- Adhesion - Adheres well to enamel, dentin.
- Bio Compatibility - mild pulp response.
- Anti Carogenic - release fluoride ions.

### Modified GIC

Metal modified GIC :- They are introduced to improve strength, fracture, toughness and resistance to wear and tear. Silver alloy is mixed.

### Resin modified GIC :-

They were introduced to overcome disturbance like moisture, sensitivity, low initial strength, working time.



## Dental Amalgam

Amalgam has been successfully used for more than a century as restorative material for tooth decay. Adv - 1.

### Classification :-

#### Based on Copper Content :-

- Low Copper alloy - less than 6% Cu
- High Copper alloy - more than 6% Cu
  - (a) Admixed
  - (b) Single Composition.

#### Based on Zinc Content :-

- Zinc containing alloy - more than 0.01% Zn
- Zinc free alloy - less than 0.01% Zn

#### Based on Shape of alloy particles :-

- Latte cut alloys
- Spherical alloys
- Spheroidal alloys.

#### Based on number of metal alloy :-

- Binary alloys - silver - tin.
- Ternary alloys - silver - tin, copper.
- Quaternary alloys - silver, tin, copper, lead.


#### Based on size of alloys :-

- Microcut.
- Macro cut.

### Setting Reaction:- (Lalke cut)

When powder is mixed with mercury, powder dissolves into it and mercury diffuses into alloy particles and plastic mass is formed. The reaction takes place in two phases:-

Stage 2 unreacted phase



### Admix:-

When powder is mixed with mercury, it diffuses into alloy particles. The mercury dissolved into Lalke cut. Particle from  $S_1$  and  $S_2$  is eliminated.



### Single Composition:-

When powder is mixed with mercury powder dissolves in mercury & mercury diffuses into powder. Silver react with mercury forming phase.



### Trituration :-

James technique is 1:1 ratio of powder. The mercury is dispersed by using a dispenser and appropriate portions of alloy is taken. The method of mixing alloy particles is known as trituration.

### Consistency :-

On trituration, Silver grey shiny mix must be obtained.

### Hand mixing :-

A mortar and pestle is used. The powder and mercury are taken in mortar and mixed using a pestle. Rotatory motion in one direction is followed.

### Mechanical mixing :-

It is done using amalgamator. Preportioned capsules are used. The capsule put on mortar. When machine is activated and holding, the capsule rotate at high speed and trituration is done. Mixing time is reduced for checking consistency. Should have shiny lustrous surface when made into balls. when made into balls when dropped should not break.

### Advantages :-

- prevent marginal leakage.
- Adequate fracture resistance.
- James long service life.

### Disadvantage :-

- Subject to corrosion and galvanic action.
- It is most aesthetic.
- They will not bind to tooth surface.

### Uses :-

- Permanent filling material in class I and II cavities.
- In retrograde root canal fillings.



## Dental Ceramics

Dental Ceramics is defined as inorganic compound with non-metallic properties typically consisting of oxygen and one or more metallic and non-metallic element that is formulated to produce the whole part of ceramic based dental prosthesis. Age - 69.

### Composition :-

- Feldspar : Basic glass former.
- Kaolin : Binders
- Alumina : Glass former and filler
- Alkaline : Glass modifiers and filler
- Colour pigment : modifies base colour.
- Opacifier : Reduce transparency.

### Classification :-

According to fusion temperature.

- High firing : 1300°C for fusion denture teeth.
- Medium firing : 1000°C for denture teeth
- Low firing : 850°C for denture teeth.
- Ultralow firing : < 850°C

### According to type :-

- Feldspar (or) Conventional porcelain
- Aluminian porcelain
- Tenite reinforced porcelain
- Glass infiltrated porcelain, Glass Ceramic.

According to use :-

- porcelain for artificial denture tooth.
- partial crown veneer
- Metal Ceramics.
- Anterior bridge Ceramics.

According to Processing method :-

- Sintering porcelain
- Cast porcelain
- Mechanical porcelain.

Advantage :-

- Highly esthetic material.
- Highly biocompatible material.
- Adequate wear resistance.
- Thermal Conductivity is very low.
- Can be modulated to any form.

Dis Advantages :-

- They are brittle and very costly.
- wear off the opposing tooth.
- very low fracture resistance.

Uses :-

- Inlays & onlays
- Single all Ceramic Crown, artificial denture
- Short Span all Ceramic bridge.

## Composite resin

### Introduction :-

It is a resin which has been strengthened by adding silica particle called fillers. Composite is a system composed of a mixture of two (or) more macromolecular essential in each other different forms.

### Uses :-

- Restoration of anterior and posterior teeth.
- To replace metal crown and fixed partial denture.
- Pit and fissure.
- Esthetic laminates.

### Type :-

Composite resins are different to classify because of the wide range of overlapping filler size.

### Based on filler part using mechanism :-

- Class I - Self cured materials.
- Class II - Light cured materials.
- Class III - Dual cured materials.

### Based on filler particle size :-

- Fine - particle size  $> 2 \mu\text{m}$ .
- ultrafine - particle size  $< 3 \mu\text{m}$ .
- microfine - Average particle size ( $0.04 \mu\text{m}$ )
- Nanofine - nano range ( $5-10 \mu\text{m}$ ).

Based on viscosity :-

- Conventional
- packable.

Restorative Composite - resin :-

Early Composite resin were generally macrofilled resin and were referred as Conventional Composite.

Supplied as :-

- Two paste System.
- powder and liquid System.

Available as :-

Single paste form in dark and light proof Syringe

Indications :-

- Resin matrix.
- filler
- Coupling Agent.





## Chemically polymerisation mechanism:

Based on mode of activation of polymerisation there are of 3 types:-

- Chemically activated resin
- Light activated resin.
- Combination Resin.

### Chemically activated resin:-

Base paste material contain benzoyl peroxide  
Catalyst paste - tertiary amine activator.

### Setting:-

when the two paste are spatulated, the amines react with the benzoyl peroxide to form free radicals.

### Disadvantage:-

- No control of curing time.
- No control of curing shrinkage.
- Hard mixing increase possibility of voids.

### Light activated Composite resins:-

Under normal light, they don't interact, however when exposed to light of correct wavelength, the photo interaction is activated.

### Dual Cure resin:

A combination of chemical and light curing is used to overcome some of the drawback. Light curing dual resin are supplied in two pastes.

### Microfilled Composites:-

The microfilled resin is introduced soon after the traditional composites. The resins achieved smoothness of unfilled acrylic filled resin.

### Nanohybrid Composite- resin:-

A new type of composite resin based on nanosized filler particles recently introduced nanohybrid.

### Bonding:-

One of the initial problem when resin restorative is microleakage which is from the shrinkage of the resin while curing.

### Acid etch technique:-

Initially developed to improve retention to enamel. At the time, it was widely practiced. Dentin could not be etched as well as enamel. Acid etching of dentin would cause injury to the pulp.

### Etchant :-

This etchant an acidic in nature. They may be grouped as

- Mineral
- Organic
- polymeric

### Specialised Composite resin :-

Flowable Composite resin.

These are hybrid Composite modified to have an increased flow.

The materials are :-

- More prone to water.
- Have half the stiffness.
- Heat polymerisation.

### Packable Composites :-

They are Composites that have very high viscosity and Surface Tension. They have higher filler loading with porous or irregular shaped particles. They are not Condensable like amalgam.

### Bioavailability :-

Internal and External toxicity of the material on the pulp & pulp involvement due to microleakage. Allergic potential on Contact to oral mucosa.

### Advantages :-

- Reliable using System allow choice of working time.
- Easy to repair
- Moderately strong and durable.
- Does not extend when compared to amalgam.


### Dis Advantages :-

- Highly technique sensitive.
- Highly shrinkable
- Stick to instrument.
- Non - Condensable like amalgam.



# RAGAS DENTAL COLLEGE & HOSPITAL

2/102, East Coast Road, Uthandi, Chennai - 600 119.



PRE CLINICAL PROSTHODONTICS RECORD NOTEBOOK

B. D. S. COURSE

NAME: J. USHA NACHIAAR

## INSTRUCTIONS

1. Students must neatly dressed and wear their coats before entering the Laboratory
2. Girls must have their hair braided and put inside the coat or have it tied-up to prevent mishap while working on Dental Motor.
3. Avoid wearing gold ornaments like rings and bangles when working with materials
4. Record books should be brought to the practical class without fail
5. Necessary instruments should be brought to the practical class by every student
6. Students should read the experimental procedure before coming to the practical class, perform the practical carefully and record the observations
7. Working place and the instruments should be kept clean
8. All instruments and equipments should be carefully handled. Any damage to the laboratory equipments must be reported to the staff in-charge
9. Materials should be used sparingly during practical work
10. The lids of the bottles and tubes should be replaced immediately after dispensing the materials
11. Burners should be closed, and motors put off when not needed
12. At the end of the every exercise, the signature of the staff member should be obtained
13. Record books should be preserved and submitted to the head of the department at the end of the academic session
14. All instruments can be procured through college stores.

# RAGAS DENTAL COLLEGE & HOSPITAL

2/102, East Coast Road, Uthandi, Chennai - 600 119.

## DEPT. OF PROSTHODONTICS

B.D.S. COURSE

PRE-CLINICAL EXERCISES

Examination Regn. No. ...54/9/1598.....

CERTIFIED that this is the bonafide record of work done by

Mr./ Miss. ....J. USHA NACHYAR.....

in the PRE CLINICAL PROSTHODONTICS, AT RAGAS DENTAL COLLEGE &

HOSPITAL, CHENNAI during the year 2021- 2022

*[Handwritten Signature]*

Signature of the Staff

*[Handwritten Signature]*  
10/6/22

*[Handwritten Signature]*  
10/6/22  
Internal

*[Handwritten Signature]*


Signature of the  
Head of the Department

## INSTRUMENTS LIST


Rubber Bowl	2 Nos
Wax knife	2 Nos
Wax carver	3 Nos
Wax spatula	3 Nos
Alginate spatula - curved spatula	1 No
Cement spatula	1 No
Upper & Lower edentulous trays (U1, U2, U3, U4, L1, L2, L3, L4)	1 (each)
Upper & Lower Dentulous trays (U1, U2, U3, U4, L1, L2, L3, L4)	1 (each)
Bonem burner with tube	1 No
Spirit lamp	1 No
Alcohol torch	1 No
Flask with Clamp	2 Nos
IC Flask	1 No
Half round file, Rat-tail file, triangular file	1 (each)
3 point articulator	2 Nos
Scissors - straight & curved	1 (each)
Camlin paint brush	2 Nos
Carbon copying pencil	1 No
Stainless steel scale	1 No
Universal plier	1 No
Glass tile (8 cm x 8 cm)	1 No
Enamel tray (square)	1 No
Small towel	1 No
Mackintosh sheet	3 Nos
Hot plate	1 No
Acrylic trimmers	4 Nos
Metal trimmers	4 Nos
PKT instruments	1 set
Plastic Knife	1 No
Straight Spatula	1 No
Porcelain mixing jar	1 No
BP handle No. 3 & No. 4	1 No
Wire-Cutter	1 No
Ceramic tile - chip syringe	1 No



# CONTENTS



EX. No.	DATE	PROCEDURE	PAGE NO.	REMARKS
02	21/11/21	PREPARATION OF VENTROCORONAL MODEL (BENDING AND REWIND)	21	
03	22/11/21	ADJUSTING LINGUAL ANCHOR IN COMPLETE MODEL	23	
04	23/11/21	EVALUATION OF CUSTOM TRAY OF EDENTULOUS MODEL	24	
05	6/12/21	Fabrication of temporary RETARD BASE (OPTIMISED PLATE) NOTES	45	
7	6/12/21	Fabrication of temporary RETARD BASE (SHELL BASE PLATE)	49	
8	11/12/21	COMPLETE DENTURE (TEETH ALIGNMENT IN CLASS 2 ARTICULATION)	53	
9	12/12/21	PROCESSING OF COMPLETE DENTURE (CURING)	69	
10	20/12/21	FINAL WORK OF COMPLETE DENTURE	77	
11				
12				

# CONTENTS

EX. No.	DATE	PROCEDURE	PAGE NO.	REMARKS
12	02/10/21	DERIVATION OF COMPLETE DENTURE (PART 1 METHOD)	25	
13	30/10/21	DERIVATION OF CAST FOR LENTHOLD'S CONDITIONED OF PARTIAL DENTURES	27	
14	05/11/21	FABRICATION OF REMOVABLE PARTIAL DENTURE (RPD)	28	



**PREPARATION OF EDENTULOUS MODELS  
(BEADING & BOXING METHOD)**

Date	Procedure	Assessment	Signature
10/10/21	Impression Maxillary Mandibular		
05/10/21	Beading of Impression Maxillary Mandibular		
05/10/21	Boxing of Impression Maxillary Mandibular		
05/10/21	Finished and polished cast Maxillary Mandibular		

# PREPARATION OF EDENTULOUS CAST

AIM:-

To prepare maxillary and mandibular edentulous cast.

APPARATUS REQUIRED:-

1. Impression compound
2. Edentulous tray
3. Type I dental plaster / Type II dental plaster
4. Modelling wax
5. Wax knife
6. Gown
7. Edentulous model
8. Rubber bowl
9. Spatula.

PROCEDURE:-

**Impression** :- A negative likeness or copy in reverse of the surface of an object, an impression of the tooth and adjacent structures of use in dentistry.

Impression compound is softened in hot water and is placed in suitable impression tray.

The edentulous model is lubricated and is pushed over the softened impression compound.

When the impression compound sets, the edentulous model is separated. Thus, the impression is obtained.

### Beading:-

When strips are completely around the periphery of both the maxillary and mandibular impressions.

The wax should be 1mm wide and 3-4 mm below the border of impression.

The logues space of mandibular impression should be also sealed with wax.

### Significance:-

The purpose of the beading is to provide the worker with a lead area which should outline and preserve peripheral wall.

### Beading:-

Beading an impression is building up vertical wall (lead) usually in wax. Producing the desired size and form of base of the cast and to preserve certain landmarks of impression.

### Preparation of alactulose cast:-


Mixed with proper H/P ratio and poured over the prepared impression.

Tapping is done to avoid air bubbles.

### Finishing and polishing:-

The cast is removed from the impression and outer surface of cast is polished with sand paper.

# ANATOMICAL LAND MARKS ON EDENTULOUS MODEL

Date	Procedure	Assessment	Signature
22/10/24	Outlining on maxillary edentulous model  Submission	5	
23/10/24	Outlining on Mandibular edentulous model  Submission	B	

## ANATOMICAL LANDMARKS

A recognizable anatomic structure used as a point of reference.

### MAXILLA

Denture structures :-

1. Labial fossae
2. Labial sulcus
3. Buccal fossae
4. Buccal sulcus
5. Hamular notch
6. Posterior palatal seal

### SUPPORTING STRUCTURES -

Primary stress bearing area:

1. Slope of residual alveolar ridge
2. postcondylar slope of hard palate

Secondary stress bearing area:-

1. Crest of residual alveolar ridge
2. Ridge

Relief area:-

1. Incisive papilla
2. Mid palatal raphe
3. Fovea palatina
4. Cephalic eminence



labial frenum:-

It is a fold of mucous membrane that extends from the labial of residual ridge to the lip.

labial sulcus:-

It occupies a potential space bounded to the labial aspect of the residual ridge. The neurolabial alveolar fold and orbicularis oris muscle. It terminates the buccal frenum.

Buccal frenum:-

It separates labial and buccal vestibule, creates the dorsal angle of the oris muscle. It is bounded by anterior orbicularis oris.

palmar: Buccinator

superior: dorsal angle of oris.

Mandibular tuberosity:-

It is a bulbous extension of residual ridge in second and third molar region. Sometimes it is double noted.

Periculae notch:-

It is a depression present between the maxillary tuberosity and ~~mandibular~~ hamulus.

Lower papilla:-

It is the exit point of maxillary nerve and vessels. It may need to be relieved to avoid pressure on nerve and vessels.

## Posterior palatal seal:-

Soft tissue along the junction of hard and soft palate which forms within physiological limits can be applied by complete denture prosthesis to aid in extension of the denture. It is identified as area between anterior and posterior seal.

functions:-

i) Retention of maxillary denture

Prevents food accumulation

Reduces gag reflex

guides position of myofascial bag.

## Residual Alveolar Ridge:-

position of alveolar process and its soft tissue covering

that remains after removal of the tooth.

## fovea palatina:-

Two small pits in the posterior aspect of palate on each side of midline or more near the attachment of soft palate to hard palate.

## Soft palatine Tapes:-

It extends from the incisive papilla to distal ends of hard palate covered with mucous membrane and has little submucosal tissue.

## Ridge:

The irregular fibrous connective ridge located in the anterior floor of hard palate.

## Mandible:

1) Dorsal structures:

1. Labial frenum
2. Labial sulcus
3. Buccal frenum
4. Buccal sulcus
5. Retromolar pad
6. Lingual frenum
7. Sublingual sulcus.

2) Supporting structures:

- primary stress bearing area:  
Buccal shelf area
- secondary stress bearing area:  
~~labial and lingual~~

## Surface Area:

1. Crest of alveolar ridge
2. Mylohyoid ridge
3. Mental foramen
4. Facial tubercle
5. Torus mandibularis.

Labial frenum:-

Contains a band of fibrous tissue that helps to attach orbicular oris muscle.

Retro-malar pad:-

Triangular pad of tissues at distal end of ridge

Buccal shelf area:-

lies at right angle to vertical occlusal force and lower teeth bounded by

medially, crest of ridge

laterally, External oblique line

Anterior, Buccal frenum

Posterior, Masseter, Retro-malar pad.

Slope of Residual Alveolar ridge

Provides more support due to underlying cancellous bone tissue to prevent a sharp, slinging. It needs to be reduced.

Mylohyoid ridge

Mylohyoid muscle attaches to the ridge, distal flange should extend below the ridge

Retro-mylohyoid cushion

posterior extension of distal portion of distal flange bounded by anteriorly, distal teeth

laterally, Pterygomandibular raphe

posteriorly, Palatoglossus surface of muscle

posteriorly, Superior constrictor muscle

Inferiorly, mylohyoid muscle and submandibular gland.

## MAXILLA



- Labial Frenum
- Labial sulcus
- Buccal frenum
- Buccal sulcus
- Residual alveolar ridge
- Maxillary tuberosity
- Fovea palatini
- Hard palate raphe
- Rugae
- Incurved papilla
- Hamular notch
- Probasal palatal seal


## MANDIBLE



Palate

- Labial frenum
- Labial sulcus
- Buccal frenum
- Buccal sulcus
- Pterygomandibular raphe
- Alveolar lingual sulcus
- Lingual frenum
- Residual alveolar ridge
- Retromolar pad
- Buccal shelf area

# FABRICATION OF CUSTOM TRAY ON EDENTULOUS MODELS (SHELLAC BASE PLATE)

Date	Procedure	Assessment	Signature
25/11/21	Marking spacer and tray outline on the cast Maxillary Mandibular	A	
25/11/21	Spacer Adaptation Maxillary Mandibular		
25/11/21	Finished custom tray Maxillary Mandibular		

# FABRICATION OF CUSTOM TRAY (Fejza adapted slough method)

## DEFINITION:-

also called as special tray, an individualized impression tray made from a cast, reinforced from a packing impression

## MATERIALS REQUIRED:-

plaster cast, wax sheet, cold mould seal, soft tin acrylic resin, acrylic hardware, pressure for polymerization

## PROCEDURE:-

- \* Make an outline on cast with a pencil
- \* In maxilla give the posterior border determined by a line extending between the buccular notch with midpoint approx. down distal to force palatine
- \* In mandible it should be just beyond the submandibular part
- \* In both by arches the outline in pressure areas should be elevated by 1mm to provide adequate space for relief landmarks
- \* Block the spare undercut with wax and adapt a layer of baseplate wax to the cast for relief of undercuts 1-2mm from the cast to the desired outline
- \* Paint the tin film on the stone cast over the relief area.

Use of the first substitute:

Removal of thin square of relief area will expose the cast thereby providing the tissue slip over the granular and Iustus region on either side.

Thin square of thin overlay

Paper-backed: The impression key material acc. to the manufacturer's recommendations. Wax and paper layers should contain checks the consistency of the resin periodically and remove when it reaches the dough stage.

Set the resin to desired thickness with a roller or a special press to make expansion material uniformly thick.

Adapt the material to cast carefully to avoid over thinning the resin on concave portions of cast.

Remove the excess key material from cast from the excessive material into handle and adapt them to key.

In maxillary key, one handle in anterior portion of the key.

In mandibular key - 3 handles as figure cast

2 on either of Iustus region

1 on anterior portion of key.

Tray thickness - 3.4mm

Height - 8mm

Height - 8mm.



place horizontal grooves across the facial and lingual  
surface of handle to improve grip

put more acrylic monomer in the expansion bag with  
collar as an eye degree at base of the attachment to  
improve bonding of the handle to bag.


Curve the expansion bag on bench or downward and  
circled plastic ball.

After setting, remove the expansion bag from the lat  
and trim the border.

Smooth the rough areas and store the bag at the last.

Ex. No. 7

## FABRICATION OF CUSTOM TRAY ON EDENTULOUS MODELS (SPRINKLE-ON METHOD)

Date	Procedure	Assessment	Signature
21/11/20	Marking spacer and tray outline on the cast Maxillary Mandibular		
21/11/20	Spacer Adaptation Maxillary Mandibular		
21/11/20	Custom tray before trimming Maxillary Mandibular	B	
21/11/20	Finished and polished custom tray Maxillary Mandibular		

# FABRICATION OF CUSTOM TRAY (Capitation method)

Definition :-

It is an individual impression tray made from a cast removed from a preliminary impression. It is used in making a final impression.

Materials required :-

primary cast, wax sheet cold mould leaf, Self cure acrylic resin, acrylic trimmer

Procedure -

- 1) Making a sulcus die, tray extension die and space die.
- 2) Sulcus die - Use a blank pencil and trace the sulcus die along the maximum depth of sulcus in preliminary cast.
- 3) Tray extension die - Use the sulcus in such a way that it follows the contour of the sulcus die but should be such that

4) Space die :-

Depend upon the type of space used, suggested spaces to used for.

- i) Maxillary arch - primary and secondary space
- ii) Mandibular arch - Modified border with 2 stops in anterior zone

2. Adaptation of wax space :-

Make sure an outline of impression tray on cast with pencil, block out the areas with wax and adapt a layer of base plate to cast for relief, wax space is given above sulcus depth.

### 2) Tissue stage:-

Wax space on either side of the ridge to a shape of either saw square or square rectangle. Cold mould seal is applied as a separating media. It increases the availability and result in even removal of wax from the die.

### 3. Deglazing:-

Light pressure on cast and wax space, which with dequid powder and dipiel is applied until there is uniform degree of approximately 2mm thick.

### 4. Fabrication of handle:-



Mix more resin, where this is dough stage, handle are formed and adapted to the Depression. Make the handle approx 3-4mm thick & 2mm height and saw height placement of horizontal groove across the facial and lingual surface of handle can improve the grip.

### 5. Trimming and Polishing:-

After setting, the die is removed from the cast and trimmed with angle of stonemason block. The areas which are rough are capable of cutting & compact which are adjusted, retained and polished.

The border of the die should be trimmed in such that it is about 1mm thick and should have smooth border without irregularities.

## FABRICATION OF TEMPORARY RECORD BASE (SHELLAC BASE PLATE)

Date	Procedure	Assessment	Signature
10/10/20	Marking record base outline on the cast  Maxillary  Mandibular	✓	
10/10/20	Finished record base  Maxillary  Mandibular	✓	

# FABRICATION OF TEMPORARY RECORD BASE (Shellac base plate)

Aim: Fabrication of temporary record base (Shellac base plate)

1. Cast
2. Bursae
3. Shellac base plate
4. Scissors
5. Collet
6. Water

Paper form of shellac base plate is selected for maxilla and mandibular arch.

Cast is soaked with slurry water before adapting the base plate. This will prevent shellac base plate from sticking to the cast while heating it.

Put the base plate over the flame, until edges place it over the cast and flame it with the bursae bursae.

Continuous finger pressure applied with wet collets adapted over the cast.

Edges the base plate and adapt it over the cast of the ridge while the material is warm and cut out the excess with scissors being applied over beyond denture areas.

Continuously heat and adapt the material to the cast.  
Using a spatula fold the material up to the surface of the  
to form a smooth rounded border.

Overheating causes shells to particles. They are porous to  
cast and stick to it.

Use files to rotate and round off the borders.

Smooth finished wax plate to be obtained as they present

Advantage:

Minimum amount of time is required for adaptation



Disadvantage:-

Chance of doing their initial adaptation since it's  
thermoelastic material

It's brittle as there is chance of breaking.

P. 110

## COMPLETE DENTURE TEETH ARRANGEMENT IN CLASS-I ARTICULATION

Date	Procedure	Assessment	Signature
21/01/21	Construction of occlusal rims  Maxillary  Mandibular	2 B	
21/01/21	Mounting on articulator		
21/01/21	Articulated casts		
21/01/21	Anterior teeth arrangement		
21/01/21	Posterior teeth arrangement		
21/01/21	Waxup & carving		



## COMPLETE DENTURE

Definition :-

A removal dental prosthesis that replaces the entire dentition and associated structures of maxilla and mandible.

Parts :-

1. Denture flange
2. Denture base
3. Denture border
4. Denture flange

Surface :-

1. Distal surface
2. Depression surface
3. Occlusal and lingual surface

Instrumentarium required :-

Impression tray, Impression compound, plaster of Paris modelling wax, self cure acrylic resin, gypsum bonded dental stone, heat cure acrylic resin, wax knife, Carver, spatula, acrylic flange, articulator, flask, clamp, bumper, rubber bowl, straight spatula, acrylic trimmer, sandpaper, cotton.

# STEPS IN COMPLETE DENTURE

## CLINICAL STEPS

1. Diagnosis and treatment planning

2. Selection of key

Primary impression using

i) Impression compound

ii) Plaster of Paris

3. Border moulding using green stick compound

4. Secondary impression with

i) Zinc oxid. eugenol

ii) Stentor's impression material

## LABORATORY STEPS

4. Beading and boxing

5. Pouring of primary cast

6. Space adaptation

7. Special key fabrication

10. Beading and boxing

11. Boxing the master cast

12. Anti-chillac wax base

13. occlusal rim fabrication

14. Jaw relation

17. Selection of felt

20. Wax strip



Pills

15. Moustache

16. Articulation

18. Teelt arrangement

19. Wax up and Carriage

21. Processing

a) flackup

b) Downup

c) Packup

d) Curup

e) Deflup

f) Lab movement

22. finishup

a) Finishup

b) Polishup

# PROCESSING OF COMPLETE DENTURE (CLASS-I)

Date	Procedure	Assessment	Signature
10/11/2024	Flasking		
	Maxillary		
	Mandibular		
10/11/2024	Dewaxing		
	Maxillary		
	Mandibular		
10/11/2024	Packing and curing		
	Maxillary		
	Mandibular		
10/11/2024	Deflasked Denture	B	<u>          </u>
	Maxillary		
	Mandibular		
10/11/2024	Finishing and Polishing		
	Maxillary		
	Mandibular		

# PRINCIPLES OF TEETH SETTING

## MAXILLA

### 1) Central incisors

front view: long axis parallel to vertical axis

side view: slightly labial inclination

occlusal view: Incisal edge is in contact with occlusal plane

### 2) Lateral incisors

front view: long axis of the tooth slightly towards the midline of the mouth

side view: Inclined labial to a greater degree than central

occlusal view: Incisal edge is same short to occlusal plane

### 3) Canine

front and side view:

long axis of the tooth is parallel to the vertical axis  
The buccal cusp half of the tooth provide its position

occlusal view: The cusp tip is in contact with the horizontal plane.

### 4) First premolar

front and side view: long axis is parallel

occlusal view: palatal cusp is same short of plane. Buccal cusp is contact with the plane.

### 5) Second premolar

front and side view: long axis parallel to the vertical

occlusal view: Both buccal and palatal cusp are in contact with occlusal plane.

vii) first molar :-

front view :- long axis slopes buccally

side view :- long axis slopes distally

occlusal view :- only the mesio palatal cusp is in contact with the occlusal plane.

viii) second molar :-

front view :- long axis slopes buccally more steeply than the first molar

side view :- long axis slopes distally more steeply than the first molar.

occlusal view :- All the cusps are chief of plane except mesio palatal cusp which is mesial to the plane.

MANDIBLE :-

i) central incisor :-

front view :- long axis parallel to vertical axis

side view :- slight labial inclination

occlusal view :- Incisal edge is seen above the plane.

ii) lateral incisor :-

front view :- long axis parallel to the vertical axis

side view :- slight labial inclination

occlusal view :- Incisal edge is seen above the plane

### ii) Canine

front view: long axis slopes very slightly towards the midline

side view: slopes slightly lingually

occlusal view: Cusp tip is 2mm above the occlusal plane

ii) first premolar:

front and side view: long axis is parallel to the vertical axis of the tooth

occlusal view: Lingual cusp is below the horizontal plane and buccal cusp is 2mm above the occlusal plane.

iii) first molar:

front view: long axis slopes lingually

side view: long axis slopes mesially

occlusally: - All the cusps are above the occlusal plane that second premolar, buccal and distal cusp are higher than mesial and lingual cusp.

iv) Second molar:

front view: Lingual inclination is more pronounced

side view: Mesial inclination is more pronounced



occlusal view: - All cusps are above the occlusal plane.

v) Second premolar:

front and side view: long axis is parallel to vertical axis


occlusal view: Both lingual and buccal cusps are 2mm above the occlusal plane.

# Relation of upper anterior teeth to Horizontal and vertical axis






Relationship of upper posterior teeth to horizontal and vertical axis



Relationship of lower anterior teeth to horizontal and vertical axis



P. 116

## REPAIR WORK OF COMPLETE DENTURE

Date	Procedure	Assessment	Signature
23/11/21	Assembling fractured parts Maxillary Mandibular		
23/11/21	Preparation of working model Maxillary Mandibular	B	<i>[Signature]</i>
23/11/21	Creation of space at fractured site and placement of retention grooves Maxillary Mandibular		
24/11/21	Processing Maxillary Mandibular		
24/11/21	Finishing and Polishing Maxillary Mandibular		

## PROCESSING OF COMPLETE DENTURE

Steps :-

- Flasking
- Removal
- Packing
- Curing
- Refinishing
- Finishing

Flasking -

The process of investing the cast and a wax replica of the desired form in a flask preparatory to mold the restorative material in to the desired products

- i) Base
- ii) lower part / Body
- iii) Lip

check the seat of trial denture to the cast and fill in deficient area with baseplate. - now take care of completely fill the guiding bands. However do not overflow the wax on the cast bands.

Select the flask that fit together smoothly without rocking place the denture and cast in the flask to check the height of denture seat in the flask

Place the dental plaster mix in the base of the flask and  
with the wax slabs and cast into wax.

Smooth the dental plaster around the cast and allow the  
stone to complete initial set.

Paint all the stone surface in base of the flask with a  
separating medium.

Place the body / Counter of flask into the position on the base  
of the flask.

Now stones into the flask, allowing line for the stone to  
flow over the denture.

Remove the stone with a finger to expose the occlusal  
surface of teeth.

A mix of dental stone is placed over the surface of the  
teeth in brief denture which is referred to a key.

Now the dental plaster over the stone surface to fill the flask.

Place the lid on the fitted flask and allow the stone to  
set.

Wax elimination

After the stone has set; The flask is placed in the  
boiling water to soften the wax for 4-6 min.

Then it is removed from water and poured, then the wax is  
washed away with boiling water.

Packing :-

Application of separating medium in the mould cavity and the cast.

For the heat cure resin is a clear mixing jar with a glass plate it inside until the resin reaches the proper stage for packing.

It is important that the mixing jar be air tight to prevent evaporation of any by products, which will cause the resin to be gummy.

After the resin has reached the dough stage, remove it from the jar, roll it and adapt to the flask.

Allow it in the bench cure for 1 hour.

Curing :-

After bench curing, place the denture in water bath at room temperature and then gradually temperature increased to 100°C for about an hour this provides polymerization better.

Short Curing Cycle :-

Processing the resin at 70°C for approx. 4 hrs and increase the temperature of the water bath to 100°C and processing for 1 hr.

Long Curing Cycle :-

Processing the denture base resin in a constant temperature water bath at 74°C for 2 hrs.

## Reflasking:

After the deaeration has been carried, the removed from the Ampy unit and allowed to bench cool.

Reflasking is the process of removal of the processed deaerated from the flask.

## Finishing, finishing and polishing:

Take the flask from the Complete Deaerated with charge bar, mounted on a laboratory bath.

finishing is done sequentially with rough and smooth grit sand paper.

Once a smooth surface is achieved polishing the deaerated with fine flume pumice in hot buff is carried out.  
data by buffing is done

Note That no pumice is employed in the intermediate and margin areas. This can be removed by using a jet of steam or bank.

pdfs

## DENTURE REPAIR

Accurate assembly and alignment of fractured denture parts in their original position  
fracture can occur in denture base and in teeth

Reason for fracture -

Denture base fracture occur due to

Alveolar resor

Harsh use

Inadequate tissue relief

Inaccurate impression

Tooth / teeth fracture occur due to

Cupped Intuferences

Faulty tooth

Detachment of retention

Constriction of angle teeth

Excessive grinding of teeth

Repairing fractured Denture:

1) Non-separated / superficial fracture:

often occur at midline frequently appearing at a deep fracture notch.

Step 2:

Broken pieces are placed together to determine accurate fit.

The denture is held together with sticky wax and small number sticks.

A latex mold of plaster is provided next. The denture after blocking the undercut.

The denture is removed from cast, free of sticky wax and sticks.

No. 558 bur to grind out fracture lines from the beginning to the end and across to increase the bonding surface. In the palate of maxillary dentures, place down-buff to strengthen the repair joint.

Coat is painted with separating medium. Place and paint autopolymerizing resin into grooves taking care to avoid air entrapment.

Build up resin repair slightly above the denture base. It is 20 min in pressure cabinet and denture fractured into two or more parts (components).

Reason:

Denture when dropped on hard surface, repair is done as similar to that done in millers fracture.

Additionally, resin cement can be used for strength.



2. fractured denture tooth with section missing

When part of denture is missing on the maxillary arch, the result of patient and lab's work, possibly repair done and inspection later. Fractured denture can be repaired on the lab and replace missing tooth.

3. Repairing denture with fractured tooth:

Method to replace the tooth depend on tooth material. It involves removal of fractured tooth and replacement.

Note:-

Align parts accurately

Do not repair III - fitting denture


Use wax - contaminated - monomer - denture repair resin

Cure denture in warm water under pressure to avoid porosity

Use wax free furnace to add rigidity.

P. 11

### REBASING OF COMPLETE DENTURE (FLASK METHOD)

Date	Procedure	Assessment	Signature
12/12/21	Preparation of working models		
	Maxillary		
	Mandibular		
12/12/21	Flasking with index		
	Maxillary		
	Mandibular		
12/12/21	Preparation of Complete dentures for rebasing		
	Maxillary	B	
	Mandibular		
12/12/21	Wax-up of Complete dentures		
	Maxillary		
	Mandibular		
12/12/21	Processing		
	Maxillary		
	Mandibular		
12/12/21	Finishing and Polishing of rebased dentures		
	Maxillary		
	Mandibular		

# REPAIRING OF COMPLETE DENTURE (fasting method)

Definition:-

The laboratory process of replacing the entire denture base material on an existing prosthesis

Indications of Repairing

Immediate denture at 3-6 months after their original construction

When residual alveolar ridges have resorbed and adaptation of denture base to the ridge is poor.

When the patient cannot afford the cost of having new denture constructed.

When the construction of new denture with the accompanying series of appointments can cause physical and mental strain

After an anterior repair of

Existing denture

Replacement material

Dentist flask

Face bow

Denture base material

clinical steps

1) static method:-

- open mouth technique (Boucher)

- closed mouth technique

2) functional method:-






chair side method

## Laboratory Diagram :-

- i) Laboratory procedures for rebaring includes artificial method, jig method and flask method.  
The method use is flask method.
- ii) The poured impression along with the denture is placed in the base flask.
- iii) A silicon mould material like is packed over the denture to involve the body.
- iv) The investment body is done for a flexible mould and flasking is complete.
- v) The flask is opened and denture is involved.
- vi) The outer denture base is trimmed leaving the first sum of angle around the procedure to retain the position of denture belt.
- vii) The belt is placed back into the mould.
- viii) The invested along perimeter the base of flask is the cut for denture.
- ix) Separating media is painted over the mould space.
- x) Plate is packed cured, finished and polished.
- xi) The finished denture are removed to check for occlusal di harmony.
- xii) Once the dentures are released, they are invested in the mould.

Ex. No. 13


# PREPARATION OF CASTS FOR KENNEDY'S CLASSIFICATION OF PARTIALLY EDENTULOUS ARCHES & SURVEYING

Date	Procedure	Assessment	Signature
21/10/21	<p>Class I</p> <p>Maxillary</p> <p>Mandibular</p> 	<p>B</p> 	
21/10/21	<p>Class IV</p> <p>Maxillary</p> <p>Mandibular</p> 	<p>B</p>	
21/10/21	<p>Surveying of the casts</p>	<p>B</p> 	


# KENNELLY'S CLASSIFICATION OF REMOVABLE PARTIAL DENTURE

## CLASS - I

Bilateral edentulous areas located posteriorly remaining natural teeth.




class of  
the lateral ectostylar area located posterior to maxillary  
nasal fossa



class II


A. Midaxial ectodermal area with neural fold  
remains both anterior and posterior to





class 2

A single but bilateral (closing the midline) edentulous  
was located anterior to the maxillary natural teeth




# FABRICATION OF REMOVABLE PARTIAL DENTURE (RPD) (INTERIM RPD)

MAXILLARY

CL - IV

MANDIBULAR

CL - I

Date	Procedure	Assessment	Signature
11/11/20	Preparation of casts interim RPD Maxillary Mandibular	B	
12/11/20	Occlusal rims Maxillary Mandibular		
13/11/20	Articulation		
14/11/20	Teeth arrangement		
15/11/20	Preparation of wrought wire clasps		
16/11/20	Wax-up		
17/11/20	Flasking and Dewaxing		
18/11/20	Packing and Curing		
19/11/20	Finishing and Polishing Maxillary Mandibular		

# REMOVABLE PARTIAL DENTURE

Definition :-

Removable partial denture defined as replacement of one or more missing teeth is a partially adentulous arch which can be readily inserted and removed from the mouth by the patient.

Instrument :-

1. Mandible and mandibular arch
2. Shellac baseplate
3. modelling wax
4. Acrylic resin
5. Articulator
6. Bunsen burner
7. Curved scissors
8. Wax spatula
9. Cold mould leaf
10. Precalcium cup
11. Acrylic resin.

Preparation of cast :-

Dentures mandible and mandibular casts are made by pouring dental stone in a mould.

Required teeth are trimmed from the cast to achieve class I or mandibular arch and class II or mandible arch.

## Fabrication of occlusal Rins:

Base plate are adapted over the edentulous ridge with adequate buccal and lingual extension or which occlusal rims are fabricated using modelling wax.

### Mounting and Articulation

The upper and lower cast with the finished occlusal rims are articulated in three point articulator according to articulation principle.

### Teeth arrangement

Artificial teeth are set, clay are made for class I mandibular cast with 19 gauge wrought stainless steel to wire.

### Wax up and Casting

Seal the base plate wax around the neck of each artificial teeth along the collar. The edentulous area are finishing & done. This is sealed with the cast. After sealing dearticulation is carried out to separate the mounted casts.

### Flasking:

It is done by reverse flasking method.

The flask should be placed in boiling water for about 10 min and then wax elimination is done.

Make sure that mould cavity is dry and a layer of separating medium is applied with a brush. Mix resin and monomer in 3:1 ratio in a clean mixing jar.

close the jar in order to prevent evaporation of moisture. When enough sludge is soaked, the material is taken out and packed.

Curing

After bench curing place the denture in boiling water bath at room temperature and then gradually the temperature is increased to  $100^{\circ}\text{C}$  for about an hour.

This provides better polymerization.

Deflasking

After the denture has been cured it is removed from the curing unit and allowed to bench cool.

Deflasking is a process of removal of pieces and dentures from the flask.

Trimming, finishing and polishing

Excess flash is removed with wheel stone trimmer. Further trimming is carried out with acrylic trimmer. Finishing is done separately with rough and smooth grade sand paper.

Once the smooth surface is achieved, polishing the denture with the fine flour pumice in wet buff is carried out later by buffing is done.

Note that pumice is not enlarged within denture. This can be removed by using a jet of clean water.

# Ragas Dental College & Hospital

No.2/102, East Coast Road, Uthandi, Chennai - 600 119.



## Pre-Clinical Record

### Conservative Dentistry

NAME SANGEETHA - G



# Ragas Dental College & Hospital

No.2/102, East Coast Road, Uthandi, Chennai - 600 119.

## DEPARTMENT OF CONSERVATIVE DENTISTRY

Examination Registration Number : 541611562

**CERTIFIED** that this is the bonafide record work done by  
**Mr./Miss** ..... SANGEETHA G .....

In the Pre-Clinical Conservative Dentistry Laboratory, at Ragas  
Dental College, Chennai.

During the year 20..17..... to 20..18.....

Signature of the Staff

Signature of  
Head of the Department

External Examiner

Date :

Internal Examiner

Date :

# Impression and Model

No.	Date	Requirements	Impression	Plaster model	Final Approval (Finishing & Polishing)	Grade
1.	16-11-17	Upper Premolar	✓	✓	✓	A-1
2.	16-11-17	Upper Premolar	✓	✓	✓	A-1
3.	16-11-17	Lower Premolar	✓	✓	✓	A-1
4.	16-11-17	Upper Molar	✓	✓	✓	A-1
5.	16-11-17	Upper Molar	✓	✓	✓	A-1
6.	16-11-17	Upper Molar	✓	✓	✓	A-1
7.	16-11-17	Lower Molar	✓	✓	✓	A-1
8.	16-11-17	Lower Molar	✓	✓	✓	A-1
9.		Lower Molar				
10.		Lower Molar				



## Exercise No. 2

### G.V. Black's Classification

Class I: ALL PIT AND FISSURE CAVITIES

- (a) Occlusal Surface of Premolars and molars
- (b) Occlusal two-thirds of facial and lingual Surface
- (c) Lingual surface of maxillary incisors.

Diagram : upper Pre Molar Class I



Line Angles : 8 LINE ANGLES

- Mesio palatal
- Disto palatal
- Mesio buccal
- Disto buccal
- Mesio pulpal
- Disto pulpal
- Bucco pulpal
- Palato pulpal

Point Angles : 4 POINT ANGLES

- Mesio buccal pulpal
- Disto buccal pulpal
- Mesio lingual pulpal
- Disto lingual pulpal

Grade :



 Signature

Diagram : Upper Molar Class I (Conservative)



Line Angles : 8 LINE ANGLES

- Mesio buccal
- Distal buccal
- Mesio palatal
- Distal palatal
- Mesio pulp
- Bucco pulp
- Linguo pulp


Point Angles : 4 POINT ANGLES

- Mesio bucco pulp
- Distal bucco pulp
- Mesio linguo pulp
- Distal linguo pulp

Grade :

  
Signature

Diagram : Upper Molar Class I (Conventional)



Line Angles : 8 LINE ANGLES

- Mesio buccal
- Distobuccal
- Mesio palatal
- Distopalatal
- Buccopulpal
- Lingopulpal
- Mesio pulpall
- Distopulpal


Point Angles : 4 POINT ANGLES

- Mesio buccopulpal
- Distobuccopulpal
- Mesiolinguopulpal
- Distolinguopulpal

Grade :

  
Signature

Diagram : Lower Pre Molar



Line Angles : 8 LINE ANGLES

- Mesio lingual
- Mesio buccal
- Mesio pulpal
- Linguo pulpal
- Distobuccal
- Distolingual
- Distopulpal
- Buccopulpal

Point Angles : 4 POINT ANGLES

- Mesio buccopulpal
- Mesio linguo pulpal
- Distobuccopulpal
- Distolinguo pulpal

Grade :

 Signature

Class II :

cavities seen on the proximal surface of posterior teeth.

Diagram : Upper Pre - Molar



Line Angles : 11 LINE ANGLES

- Faciopulpal
- Linguopulpal
- Distopulpal
- Axio pulpals
- Distofacial
- Distolingual
- Faciogingival
- Axioingival
- Axiofacial
- Axio lingual


Point Angles : 6 POINT ANGLES

- Distofaciopulpal
- Distolinguo pulpals
- Axio bucco pulpals
- Axio buccogingival
- Axio linguopulpal
- Axio linguogingival

Grade :

Signature

Diagram : Class II - Upper First Molar (MO)



Line Angles : 11 LINE ANGLES

- Distobuccal
- Distopulpal
- Distolingual
- Linguspulpal
- Axio-pulpal
- Axio lingual
- Axiofacial
- Faciogingival
- Linguogingival
- Axio-buccal
- Buccopulpal


Point Angles : 6 POINT ANGLES

- Distobuccalpulpal
- Distolinguspulpal
- Axio-buccalpulpal
- Axio-linguspulpal
- Axio-linguogingival
- Axio-buccogingival

Grade :

Signature

Diagram : Class II - Upper Lower Molar (MO)



Line Angles : 11 LINE ANGLES

- Distofacial
- Axiofacial
- Axio lingual
- Distolingual
- Axio gingival
- Facio gingival
- Linguo gingival
- Facio pupal
- Linguo pupal
- Axio pupal
- Distopupal


Point Angles : 6 POINT ANGLES

- Distofaciopupal
- Axiofaciopupal
- Axiofacio gingival
- Axio linguo gingival
- Axio linguo pupal
- Distolinguopupal

Grade :

Signature

Diagram : Class II for amalgam - Upper Lower Molar (MO)



Line Angles : 11 LINE ANGLES

- Disto facial
- Axio lingual
- Disto lingual
- Axio facial
- Faciogingival
- Axioingival
- Linguogingival
- Faciopulpal
- Linguopulpal
- Axio pulpall
- Distopulpal

Point Angles : 6 POINT ANGLES

- Disto facio pulpall
- Axio facio pulpall
- Axio facio gingival
- Axio linguogingival
- Axio inguopulpall
- Distalinguopulpall


Grade :

Signature



Class III : Cavities seen on the provisional surface of anterior teeth which do not involve incisal edges.

Diagram : Class III for Composite - Upper Central Incisor



Line Angles : 6 LINE ANGLES

- Faciogingival
- Linguogingival
- Axioingival
- Axiofacial
- Axio lingual
- Incisal

Point Angles : 3 POINT ANGLES

- Axiofaciogingival
- Axio linguogingival
- Axio incisal

Grade :

  
Signature

Class IV :

Diagram : Class IV for Composite - Lower Central Incisor



Line Angles : 11 LINE ANGLES

- Mesiofacial
- Mesiolingual
- Mesio-pulpal
- Axiofacial
- Axio-lingual
- Axio-lingival

Point Angles : 6 POINT ANGLES

- Mesiofacio-pulpal
- Mesio-linguo-pulpal
- Axiofacio-pulpal
- Axio-linguo-pulpal
- Axiofacio-lingival

Grade :

- Axio-pulpal
- Facio-lingival
- Linguo-lingival
- Linguo-pulpal
- Facio-pulpal
- Axio-linguo-lingival

  
Signature

Class V :

Diagram : Class V - Upper Molar



Line Angles : 8 LINE ANGLES

- Mesioincisal
- Mesioingival
- Distoincisal
- Distolingival
- Axio mesial
- Axio distal
- Axio incisal
- Axio lingival

Point Angles : 4 POINT ANGLES

- Axio mesio lingival
- Axio disto lingival
- Axio mesio incisal
- Axio disto lingival

Grade :

 Signature

Diagram : Class V for Lower Molar



Line Angles : 8 LINE ANGLES

- Mesioincisal
- Mesiogingival
- Distoincisal
- Distogingival
- Axio mesial
- Axiodistal
- Axioincisal
- Axio gingival

Point Angles : 4 POINT ANGLES


- Axio mesioincisal
- Axio disto incisal
- Axio mesio gingival
- Axio disto incisal

Grade :

 Signature

#### STEP 4: CONVENIENCE FORM

Shape of the Preparation that provides for adequate desolation accessibility and ease of operation in preparing and out storing the tooth.



#### STEP-5 - REMOVAL OF ANY REMAINING INFECTED DENTIN (OR) RESTORATIVE MATERIAL

Elimination of any infected carious tooth structure or faulty restoration material left in the tooth after initial tooth preparation.

#### FEATURES -

It affects the aesthetic result of new restoration.

Note retention is needed for new restoration.

Radiographic evidence of caries under old restoration.

Pulp symptomatic pre-operatively.

Peripherally (or) remaining old restorative material is not needed.

#### STEP 6: PULP PROTECTION

Pulp should be protected with either cavity varnish liners or base.

zinc phosphate cement being a high strength base with compressive strength of 104 mpa is preferred for pulp protection.

#### STEP 7: SECONDARY RESISTANCE AND RETENTION FORM

Compound and complex preparation require additional retentive, resistance features.

## Features: -


- (i) Mechanical features
  - (a) Retention locks grooves
  - (b) Preparation retention
  - (c) Skirts
  - (d) Beveled chamfer margins
  - (e) Fine and slope
- (2) placement of etchers, primer (or) adhesive on prepared wall.




## STEPS: - FINISHING OF EXTERNAL WALLS.

Development of a specific cavosurface design and degree of smoothness (or) roughness that produce the maximum effectiveness of restorative material being used.

## STEP 9: - FINAL PROCEDURES.

- The final procedure of tooth preparation involve
- (i) Cleaning and inspecting
  - (ii) Desensitizing
  - (iii) Sealing of the dentinal tubules.



Name of Instrument	Diagram	Description & Uses
Tweezer		<ul style="list-style-type: none"> <li>* Consists of two arms joined by one end and the other remains apart</li> <li>* Used for carrying cotton rolls and wedges in oral cavity</li> </ul>
Excavator		<ul style="list-style-type: none"> <li>* Used to excavate cavities</li> <li>* Used to remove undetermined enamel</li> <li>eg. Hatchet, hoe, discoid spoon etc.</li> </ul>
Enamel Chisel		<ul style="list-style-type: none"> <li>* Blade is in the same length as is parallel to it</li> <li>* Used for clearing or splitting of undetermined enamel.</li> </ul> <p>Formula 10-7-14</p>

Hoe



- \* Primary cutting edge is perpendicular to axis of handle.
- \* Used for smoothing cavity floor, planning cavity preparation, forming angles.
- \* Pull/pinch/downward cutting motion.
- Eg formula  $4\frac{1}{2} - 1\frac{1}{2}$ .

Hatchet






- \* Bilateral Instrument
- \* Cutting edge parallel to long axis of handle.
- \* Used for removing unsupported enamel, undetermined enamel.
- Formula  $3 - 2 - 2\frac{1}{2}$ .




Gingival Margin Trimmer









- \* Modified hatchet
- \* Double flanged instrument
- \* Rounding of axis
- pulpal length
- angle of two surface preparation








Name of Instrument	Diagram	Description & Uses
Glass Slab & Cement Spatula		<ul style="list-style-type: none"> <li>★ Used for mixing powder and liquid for base</li> <li>★ Used for alloys setting</li> <li>★ Used for preparing homogenous mix of liquid and powder.</li> </ul>
Plastic Spatula		<ul style="list-style-type: none"> <li>★ Used to mix cements like Gac, zinc, poly carboxylate</li> <li>★ Easy to clean</li> <li>★ Cement do not stick to spatula</li> </ul>
Plastic Instrument		<ul style="list-style-type: none"> <li>★ Carrying mixed cement</li> <li>★ Used to place base prepared cavity</li> <li>★ Used to cements which may abrad spatula</li> <li>★ Tend to stick to spatula</li> </ul>


Name of Instrument	Diagram	Description & Uses
Mortar & Pestle		<ul style="list-style-type: none"> <li>★ Mortar and pestle for amalgam and repeated grinding till it is absorbed</li> </ul>
Amalgam Carrier		<ul style="list-style-type: none"> <li>★ Used in filling the restoration</li> <li>★ Helps in carrying amalgam to prepare tooth surface</li> </ul>
Amalgam Condenser  (i) Serrated (Cylindrical)  (ii) Serrated (Parallelogram)		<ul style="list-style-type: none"> <li>★ Used to condense amalgam into casting</li> <li>★ Facilitate contouring buccal and lingual surface</li> <li>★ Prevents slippage of amalgam</li> </ul>

Name of Instrument	Diagram	Description & Uses
Amalgam Burnisher (i) Ball Burnisher (ii) Conical Burnisher (iii) T-Burnisher		<ul style="list-style-type: none"> <li>Used for smoothing, finishing, polishing of restoration</li> <li>Helps to give good appearance to the repaired tooth</li> </ul>
Amalgam Carver a) Hollenbeck b) Wards c) Frahms		<ul style="list-style-type: none"> <li>For removing excess amalgam</li> <li>For shaping amalgam to natural anatomic contour</li> <li>To show cavo surface margin</li> <li>To curve the morphology of tooth restoration</li> </ul>
Wedges		<ul style="list-style-type: none"> <li>Made of wood, metal, celluloid or plastic</li> <li>To stabilize matrix band</li> </ul>

Name of Instrument	Diagram	Description & Uses
Matrix Retainer & Band 1. Ivory No. 1		<ul style="list-style-type: none"> <li>Used to hold Matrix to provide wall for proximal surface.</li> <li>Easier to place unilateral partial coverage band.</li> </ul>
2. Ivory No. 8		<ul style="list-style-type: none"> <li>Provide band for encircling either crown of tooth surface for class II and MOD activities.</li> <li>Used for No MOD and DO.</li> </ul>
3. Tofflemire		<ul style="list-style-type: none"> <li>Universal Matrix Retainer</li> </ul>

Name of Instrument	Diagram	Description & Uses
Bur/Diamond abrasives 1.245		<ul style="list-style-type: none"> <li>* In head &amp; cemented carbide bur, particles are held together.</li> <li>* Uses :- In cutting cavity.</li> </ul>
2. Round		<ul style="list-style-type: none"> <li>* Spherical shape extension of preparation of retention features.</li> <li>* Used in class II - molar</li> </ul>
3. Straight Fissure		<ul style="list-style-type: none"> <li>* It is an elongated cylindrical bar.</li> <li>* Used for amalgam cavity preparation.</li> </ul>
4. Tapering Fissure		<ul style="list-style-type: none"> <li>* Slightly tapered cone cut to small end directed away from shank.</li> <li>* Uses : For indirect retention.</li> </ul>

Name of Instrument	Diagram	Description & Uses
Handpiece		<p>It is the instrument which transmit power to rotary cutting or abrading instrument.</p> <p><u>Classification</u></p> <ul style="list-style-type: none"> <li>→ Straight handpiece.</li> <li>→ Contra angled hand piece.</li> <li>→ Right angled hand piece.</li> </ul> <p><u>Based on head design</u></p> <ul style="list-style-type: none"> <li>→ Standard</li> <li>→ Mini</li> <li>→ Torque</li> </ul> <p><u>Maintenance:-</u></p> <ul style="list-style-type: none"> <li>→ Sterilising in autoclave.</li> <li>→ Periodic oiling of hand piece.</li> <li>→ oiling is done with help of pressurized container.</li> <li>→ Surface disinfection of contamination of handpiece.</li> </ul>

Name of Instrument	Diagram	Description & Uses
Handpiece		<p>It is the instrument which transmit power to rotary cutting or abrading instrument.</p> <p><u>Classification</u></p> <ul style="list-style-type: none"> <li>→ Straight handpiece</li> <li>→ Contra angled hand piece.</li> <li>→ Right angled hand piece.</li> </ul> <p><u>Based on head design</u></p> <ul style="list-style-type: none"> <li>→ Standard</li> <li>→ Mini</li> <li>→ Torque</li> </ul> <p><u>Maintenance:-</u></p> <ul style="list-style-type: none"> <li>→ Sterilising in autoclave</li> <li>→ Periodic oiling of hand piece.</li> </ul>

## Exercise 5 :

### Cavity Preparation in Plaster Models & Wax Up

S.No.	Date	Exercise	Tooth	Cavity Prep	Wax Up
1.	22-1-18	Class I	Upper Pre Molar	—	Check
2.	22-1-18	Class I	Lower Pre Molar	—	Check
3.	1-2-18		Upper Molar - (Conservative)	—	Check
4.			Upper Molar (Conservative)		
5.	22-3-18	Class I	Upper Molar (Palatal Extn.)	} <i>[Handwritten notes]</i>	Check
6.	22-2-18	Class I	Lower Molar (Buccal Extn.)		Check
7.	12-3-18	Class II	Upper Pre Molar (MO)	✓	Check
8.	19-2-18	Class II	Lower Molar (MO)	✓	Check
9.		Class II	Lower Molar (MOD)		
10.		Class II	Upper Molar (MO)		



## Exercise 6 : Natural Tooth

S.No.	Date	Tooth No.	Requirements	cavity	Base & matrix	Restoration
1.	16.4.18	14	Class I	PMU	Direct	
2.	23.4.18	34	Class I	PMU	Direct	
3.	19.4.18	16	Class I (Conservative)	PMU	Direct	
4.	19.4.18	26	Class I (Conventional)	PMU	Direct	
5.	15.5.18	46	Class I (Buccal extn.)	PMU	Direct	
6.	6.6.18	17	Class I Palatal Extn.	PMU	Direct	
7.	20.7.18	27	Class II (MO)	PMU	Direct	
8.	20.7.18	15	Class II (DO)	PMU	Direct	
9.	20.7.18	36	Class II (MO)	PMU	Direct	
10.	20.7.18	37	Class II (MOD)	PMU	Direct	
11.		11	Class III			
12.		21	Class IV			

## Exercise 7 : TYPHODONT

S.No.	Date	Requirements	Tooth No.	Cavity Preparation	Base & Matrix	Restoration
1.	2-7-18	Class I	14	mod	mod	
2.	2-7-18	Class I (Conservative)	16	mod	mod	
3.	21-7-18	Class I (Conventional)	17		mod	18/8/18
4.	21-7-18	Class I (Buccal Extn)	37	high 20/12	mod	
5.	21-7-18	Class I (Palatal Extn)	27		mod	
6.	21-7-18	Class I	34		mod	
7.	21-7-18	Class II (MO)	15		mod	
8.	21-7-18	Class II (DO)	26	high 20/12	mod	
9.		Class II (DO)	35		mod	
10.	21-7-18	Class II (MO)	36		mod	
11.	21-7-18	Class II (MOD)	46		mod	
12.		Class III	11			
13.		Class IV	21			
14.		Class V	13/23			

# Sterilization of hand instruments

INSTRUMENT	METHOD	STORAGE ADVANTAGE	DIS ADVANTAGE
Restorative and Endodontic Instrument	Autoclave Chemical Vapour pressure-sterilization	Sterilising Clothe. Surgical packs	<ul style="list-style-type: none"> <li>→ Item sensitive to high temperature</li> <li>→ cannot be autoclave</li> <li>→ steam can corrode steel on carbide bur</li> </ul>
Carbon Steel Instrument burs	Dry heat and chemical vapour pressure	Carbon steel and Corrode sensitive burs on instrument	<ul style="list-style-type: none"> <li>→ high temperature can damage sensitive towels and clothes</li> </ul>
Hand piece	Autoclave Chemical Vapour Steriliza tion, Bkays no oxide gas by circulating autoclave for min at 12°C	Carbon steel instrument bur do not corrode	<ul style="list-style-type: none"> <li>⇒ High temper- ature damage rubber and plastic</li> </ul>
Files and Endo- dentic instru- ments	Heat or Soft sterili- zation, Heat vapour method works like autoclave but it was alcohol.	Heat can kill blood pathogen in place when sterilized	Doesn't kill spores, can't be used for surgical instrument

## Exercise 9 :

### Material Used in Practical Exercises :

Materials	ADA Specification	Composition	Chemical Reaction	Uses
Zinc Phosphate Cement	ADA no: 8	<p>Powder: - Zinc oxide powder = 90%</p> <p>LIQUID: - Phosphoric acid = 33.2% Water = 36% Aluminium Phosphate = 12% Zinc = 7.2%</p>	<p>Acid base reaction Phosphoric acid attacks Zn ions The Aluminium complex reacts with <math>H_3PO_4</math> to form <math>H_2PO_4^-</math> gel.</p>	<p>As base temporary restoration Lining of restoration and orthodontic bands.</p>
Zinc Oxide Eugenol	ADA no: 30	<p>Zinc oxide white resin Zinc stearate, zinc acetate, <math>MgO_2</math>, Eugenol, olive oil clove oil</p>	<p><math>ZnO</math> hydrolyzed followed by reaction between <math>OH_2</math> and eugenol to form acetate</p>	<p>As lined and bases orthodontic cementation Lining permanent restoration</p>
Zinc Polycarboxylate	ADA no: 61	<p><math>ZnO</math>, <math>MgO</math> Stannous fluoride <math>H_2O</math>, polyacrylic acid, it's tartaric acid</p>	<p><math>Zn</math>, <math>Mg</math> <math>Sn</math> are released by acid</p>	<p>As base temporary filling material Temporary permanent cementation</p>

Materials	ADA Specification	Composition	Chemical Reaction	Uses
Glass Ionomer cement	ADA No. 66	<p>Resin Fluoride Aluminium Silicone</p> <p>Liquid polyacrylic acid</p>	Acid base reactions	<p>Amalgam replacement</p> <p>Restoring agent</p> <p>• pit and fissure sealant</p>
Amalgam	ADA No. 1	Alloy made by mixing mercury and Silver, tin alloy	$Ag_2Sn + Hg \rightarrow Ag_2Hg_2 + Sn_2Hg + Ag_2Sn$ (unreacted)	<p>1) For restoration of cavities</p> <p>2) permanent filling material in class I, II, III restoration all core material</p>

Dr. [Signature]  
20/6/18